Instructions for use: P. 2

POLARIS® VALVE

Adjustable pressure valve for CSF shunting

Sterile, single use

&

POLARIS® VALVE ADJUSTMENT KIT

Device for reading and adjusting the operating pressure of the Polaris® valve

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
Instructions for use

1. Indications	6
2. Contra-indications	6
3. Description and Operating principle of Polaris® valves	6
4. Description and Operating Principle of the Polaris® Adjustment Kit	7
5. Configurations of the Polaris® Adjustable Pressure Valve	8
6. Measurement Unit and Calibration of Operating Pressures	8
7. Behavior during Magnetic Resonance Imaging (MRI)	9
8. Sterilization – Decontamination of Polaris® Valves and Valve Kits	10
9. Cleaning – Decontamination of the Polaris® Adjustment Kit	10
10. Instructions	10
Choice of valve model	10
Adjustment of a Polaris® valve BEFORE IMPLANTATION: selection of operating pressure	10
Implantation Technique	12
Post-operative X-ray examination: Identification of the Valve Model and Pressure Reading	13
Patency Test (post-operative)	13
Sampling the CSF and injection	14
Adjustment of a Polaris® valve AFTER IMPLANTATION: reading and/or changing the operating pressure setting	14
Potential Causes of Adjustment Difficulties and Troubleshooting	16
11. Precautions for the Daily Life of the Patient	17
12. Complications / Side effects	17
Obstruction	17
Infection	17
Overdrainage	17
Other	18
13. Guarantee	18
14. Processing of the products after use	18
Destruction after use	18
Return of products	18
15. Symbols	18
References	19

Figure 1:

- Polaris® adjustable pressure valve (SPV, SPVA, and SPVB models). Top and side views

- *Polaris® adjustable pressure valve (SPV-SX, SPVA-SX and SPVB-SX models). Top and side views*



Figure 2: Composition of a complete Polaris® kit (SPV-2010 model)

WARNINGS:

FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

READ THE INSTRUCTIONS FOR USE CAREFULLY BEFORE ADJUSTING THE VALVE OPERATING PRESSURE.

1. Indications

The Polaris® adjustable pressure valve is designed for the treatment of hydrocephalus by shunting the Cerebrospinal Fluid (CSF) to the abdominal cavity or right atrium of the heart.

The Polaris® valve adjustment kit is designed for reading and adjusting the operating pressure of a Polaris® valve in order to adapt to the clinical evolution of the patient.

2. Contra-indications

The contra-indications are the following:

- established or suspected infections along the length of the shunt (meningitis, ventriculitis, peritonitis, septicemia or bacteremia) or any infection present in any part whatsoever of the body
- patients on anticoagulant therapy, or presenting with bleeding diathesis, or with hemorrhagic CSF (indeed, the presence of blood in the system could lead to an obstruction in the system)
- ventriculo-atrial shunts in patients with congenital cardiopathies or other malformations of the cardio-pulmonary system


However, if the clinical benefit outweighs these contra-indications, the implantation of a shunt in such cases is performed under the responsibility of an experienced neurosurgeon. The clinical condition of the patient must, therefore, be subject to increased surveillance.

WARNING:

DO NOT USE AN EXTERNAL SHUNTING DEVICE (DRAINAGE BAG, ETC.) IN SERIES WITH A VALVE AS THE TWO SYSTEMS COULD INTERFERE WITH EACH OTHER AND DISRUPT CONTROL OF THE DRAINAGE.

3. Description and Operating principle of Polaris® valves

Figure 3: The Polaris® adjustable pressure valve

A precision valve

The Polaris® adjustable pressure valve (Figures 1&3) is a single use implantable device.

The Polaris® valve allows drainage of Cerebrospinal Fluid (CSF) in a given direction.

The CSF arrives in the valve through the inlet connector [1], passes through the body of the valve [5] and leaves through the outlet connector [6].

The connectors are made of stainless steel and the body of the valve is in polysulfone.

On either side of the inlet connector [1] are suture holes [2] which make it possible to attach the valve to the sub-cutaneous tissues to prevent the valve from migrating.

On the upper surface of the valve an arrow [4] shows the direction of CSF flow through the valve. This helps to position the valve correctly during implantation.

On the lower surface of the valve there is an individual serial number.

The valve body contains a ball-in-cone mechanism which determines the operating pressure of the valve.

In conditions of normal use this mechanism provides an anti-reflux function and is not sensitive to temperature variations.

The valve body, which cannot be deformed, protects the mechanism from shocks.

It also prevents attempts at pumping or puncturing the valve and makes it insensitive to variations in percutaneous pressure.

The Polaris® valve is a latex-free and phthalate-free product.

Adjustable pressure

The resistance of the Polaris® valve can be adjusted without re-operating, in order to adapt it to the clinical evolution of the patient.

The principle is based on the pressure variation exerted on a ruby ball [7] by a flat, semi-circular spring [8] at different points on its curvature.

The spring is attached to a rotor [9] which is able to rotate inside the valve body on its ruby central axis [11].

The operating pressure of the Polaris® valve is determined by the angular position of its rotor.

There are 5 pressures available. These correspond to the 5 possible positions of the rotor.

Position No. 1 corresponds to the lowest pressure and Position No. 5 to the highest pressure.

Titanium radio-opaque studs [10] inserted on the right-hand side of the valve body indicate the 5 positions of the rotor/spring assembly (cf. §10 - Post-operative X-ray examination).

The Polaris® adjustable pressure valve exists in 4 different pressure ranges (cf. Table of Operating Pressures).

To the left of the inlet connector there are radio-opaque points [3] which identify the pressure range of the valve (Figures 3&16, and §10 - Post-operative X-ray examination).

Table of Operating Pressures for Polaris® valves

References concerned	SPV-140 SPVA-140	SPV, SPV-2010 SPVA, SPVA-2010 SPVB, SPVB-2010	SPV-300 SPVA-300	SPV-400 SPVA-400
Identification of the range	0 point	1 point	2 points	3 points
Available pressures (mmH ₂ O)	1 10	30	50	80
	2 40	70	100	150
	3 80	110	150	230
	4 110	150	220	330
	5 140	200	300	400

Magnetic lock

The rotor of the Polaris® valve is fitted with a patented magnetic locking system.

This self-locking system is based on the permanent reciprocal attraction of 2 mobile micro-magnets of opposite polarity.

This “magnetic lock” holds the rotor in the selected position, thus preventing any accidental change in operating pressure if the valve is exposed to magnetic fields.

In fact, in the presence of a standard magnetic field (unidirectional) the two micro-magnets are attracted in the same direction, so only one of the two magnets moves in the direction of the field, while the other remains locked.

Changing the operating pressure of the valve first requires the simultaneous unlocking of the two micro-magnets in the valve by a specific magnet.

The rotor can then turn freely on its central axis.

4. Description and Operating Principle of the Polaris® Adjustment Kit

Changing the operating pressure of the valve is done percutaneously using an adjustment kit specifically designed to adjust the Polaris® valve. It has 3 parts:

Locating Instrument (Locator)

Reading Instrument (Compass)

Setting Instrument (Magnet).

(cf. Table of References and Figures 4, 4a, 4b, and 4c):

Figure 4: The Polaris® adjustment kit

The components of the Polaris® adjustment kit are packaged, non-sterile, in a re-usable box (Figure 4).

The Polaris® adjustment kit is a latex-free product.


Locating Instrument

The Locating Instrument (**Locator**, Figure 4a) is designed to locate the valve under the skin, an essential step in reading and adjusting the operating pressure.

The Locator houses the Reading Instrument and the Setting Instrument in turn.

It enables Polaris® valves to be adjusted, whatever their pressure range, from among the 4 ranges available (cf. §3 - Table of Operating Pressures).

Figure 4a: Locating Instrument (Locator with three special pressure reading rings)

The reading area [1] on the Locator shows the reference of the SPV valve and the 5 “standard” pressures relating to this product reference.

The reading areas [2] on the reading rings show the “special pressure” valve

references, SPV-140, SPV-300 and SPV-400, as well as each of the 5 pressures relating to these references (cf. §3 - Table of Operating Pressures).

Two pins [3] on the rim of the Locator make it possible to attach the reading rings so that the Locator will display the pressure range relating to the valve model to be adjusted.

To adjust a valve, the Locator must be correctly positioned above the implant.

The green arrow [4] shows the direction of CSF flow, and makes it possible to position the Locator correctly in relation to the valve.

The Locator center cut-out area in the shape of an imprint of the valve [5] is used to position the Locator as close as possible to the implant.

Two guide pins [6] enable the Reading Instrument to be correctly positioned inside the Locator.

The cut-out in the shape of an imprint of the Setting Instrument [7] makes it possible to position the Magnet correctly inside the Locator.

Reading Instrument

The Reading Instrument (**Compass**, Figure 4b) fits inside the Locator and allows a reliable, accurate, direct reading to be taken of the valve operating pressure.

The Compass determines the position of the valve rotor. This corresponds to the operating pressure that has been read.

Figure 4b: Reading Instrument (Compass)

The operating pressure reading is determined by the alignment of the Compass needle [1] with one of the lines on the contour of the Compass [2] and the corresponding pressure value in the reading area of the Locator.

To position the Compass correctly in the Locator, the two lateral guide pins [3] of the Compass must be inserted into the guide pins on the Locator.

Setting Instrument

The Setting Instrument (**Magnet**, Figure 4c) is made up of a specific assembly of magnets designed to enable the valve rotor to be unlocked and the operating pressure of the Polaris® valve to be adjusted.

Figure 4c: Setting Instrument (Magnet)

The Magnet is inserted into the Locator.

There are two lateral lugs [1] to enable it to be positioned correctly.

The insertion of the Magnet into the Locator makes it possible:

- to unlock the valve
- to adjust the operating pressure by turning the Magnet.

Clockwise rotation increases the operating pressure of the valve, while counter-clockwise rotation reduces it.

A green marker [2] pointing to the pressure values in the Locator reading area and the “clicks” produced when the Magnet rotates make it possible to determine the successive operating pressures set by the Magnet.

Removing the Magnet activates the automatic re-locking of the valve rotor due to the mutual attraction of the two mobile micro-magnets.

PRECAUTIONS:

DO NOT STORE OR HANDLE THE MAGNET NEAR ANY DEVICE LIKELY TO BE AFFECTED BY A POWERFUL MAGNETIC FIELD.

DO NOT PLACE THE MAGNET CLOSE TO A POWERFUL MAGNETIC FIELD (E.G.: MRI). IT COULD BECOME A DANGEROUS PROJECTILE, OR ITS PERFORMANCE COULD BE ADVERSELY AFFECTED.

HANDLE THE MAGNET WITH CARE NEAR METALLIC OBJECTS (E.G. HOSPITAL FURNITURE). THERE COULD BE A VERY FORCEFUL ATTRACTION BETWEEN THEM.

5. Configurations of the Polaris® Adjustable Pressure Valve

The Polaris® adjustable pressure valve is available in 4 different pressure ranges, each of which offers 5 pressures (cf. §3 - Table of Operating Pressures).

The Polaris® valves in the “standard” pressure ranges (30-200mmH₂O) are available in 3 models:

- without a reservoir (SPV)
- with an integrated reservoir of the antechamber type (SPVA)
- with an integrated reservoir of the burr hole type (SPVB).

These 3 models are available as the valve alone or as a complete kit (Figures 1 & 2).

The Polaris® valves in the “special” pressure ranges (10-140, 50-300, and 80-400mmH₂O) are available in 2 models:

- without a reservoir (SPV-140, SPV-300 and SPV-400)
- with an integrated reservoir of the antechamber type (SPVA-140, SPVA-300 and SPVA-400)

These “special pressure” valves are not available as complete kits; only the valve is available.

The standard pressure Polaris® valves (SPV, SPVA and SPVB), as well as the special low pressure valves (SPV-140 and SPVA-140), are also available with a SiphonX® “gravitational” anti-siphon device, pre-connected downstream of the valve.

SiphonX® thus makes it possible to add additional resistance to the operating pressure of the Polaris® valve. This resistance depends on the position of the patient (cf. §6 - Table of Operating Pressures for Polaris® valves with pre-attached SiphonX®).

The 5 corresponding models (SPV-SX, SPVA-SX, SPVB-SX, SPV-140-SX and SPVA-140-SX) are only available as the valve alone form (with pre-attached anti-siphon device) and are not available as complete kits.

Sophysa offers a complete range of ventricular and distal radio-opaque catheters which allow the CSF to flow to the valve and from the valve to the peritoneum or right atrium respectively, depending upon the type of shunt chosen by the neurosurgeon. A complete Polaris® shunt system must consist of a ventricular catheter, a Polaris® valve and a distal catheter (atrial or peritoneal).

Patient Identification Card

Each Polaris® valve is packaged with a Patient Identification Card (PIC) and traceability labels.

The surgeon is responsible for filling in this Patient Identification Card and giving it back to the patient.

6. Measurement Unit and Calibration of Operating Pressures

The pressures mentioned are in mmH₂O.

1 mmH₂O corresponds to 9.807 Pa or 0.074 mmHg.

The valves are calibrated on the basis of a flow rate of 10ml/h.

Each valve is tested individually: the measurement concerns the upstream pressure of a 10ml/h flow of water passing through the valve and the Sophysa proximal and distal catheters.

The calibration is performed disregarding the resistance of the catheters.


Thus the pressures given on the valve labels correspond to the resistance of the valve alone.

The catheters add their own resistance to the shunt.

The reservoirs are not considered to have any particular resistance.

The calibration of SPV valves is performed with a tolerance of -10/+15mmH₂O on the measured pressures.

Figure 5: Flow Rate - Pressure Curves for the SPV Model

This curve is obtained by varying the applied pressure for each pressure setting and measuring the resulting flow rate. The values are given disregarding the resistance of the catheters.

Table of mean values at 10ml/h for each position setting of the Polaris® valve, and corresponding values at 20ml/h:

SPV Polaris® valve position setting	at 10ml/h (mmH ₂ O)	at 20ml/h (mmH ₂ O)
Position 1	30	35
Position 2	70	75
Position 3	110	115
Position 4	150	160
Position 5	200	225

The effects on the device operating pressure of changes in the position of the patient and of sub-cutaneous pressure are negligible.

Case of a valve preconnected to a SiphonX® gravitational anti-siphon device

In this case, the position of the patient has a direct influence on the system operating pressure (cf. Table hereafter):


Table of Operating Pressures for Polaris® valves with pre-attached SiphonX®

References concerned		SPV-140-SX SPVA-140-SX		SPV-SX, SPVA-SX, SPVB-SX	
SiphonX® orientation		0°	90°	0°	90°
					
Available pressures (mmH ₂ O)	1	10	210	30	230
	2	40	240	70	270
	3	80	280	110	310
	4	110	310	150	350
	5	140	340	200	400

When the "Valve + SiphonX®" assembly is horizontal, SiphonX® does not add any additional resistance to the operating pressure of the Polaris® valve, and the assembly therefore behaves like a valve on its own (cf. Figure 5).

When the "Valve + SiphonX®" assembly is vertical, SiphonX® adds 200mmH₂O to the operating pressure of the Polaris® valve (cf. Figure 6).

Figure 6: Flow Rate - Pressure Curves for the SPV-SX Model (SPV with a pre-attached gravitational anti-siphon device SiphonX®) in vertical position.

This curve is obtained by varying the applied pressure for each pressure setting and measuring the resulting flow rate.


The values are given disregarding the resistance of the catheters.

Thus the pressures given on the labels of SPV-SX model correspond to the resistance of the "Valve + SiphonX®" assembly.

The catheters add their own resistance to the shunt.

For all intermediate inclinations between the horizontal and the vertical, the SiphonX® adds a resistance which depends on the angle of inclination (cf. Figure 7).

Figure 7: Operating pressures for the SiphonX®, SX-200 model, depending on its inclination, whatever the flow rate between 5 and 50ml/h.

This curve is obtained by applying the following method: the measurement concerns the upstream pressure of a 10ml/h flow of water passing through the SiphonX® and the Sophysa proximal and distal catheters and by varying the angle from 0° (horizontal) to 90° (vertical). The same curves are obtained for flow rates between 5 and 50ml/h.

The measurement is performed disregarding the resistance of the catheters.

7. Behavior during Magnetic Resonance Imaging (MRI)

A shunt consisting of a Polaris® valve (including connectors, and possibly reservoirs and SiphonX® gravitational anti-siphon device) and its catheters, is considered as "MR Safe" in accordance with the definition in the standard, ASTM F-2503-05.

A patient fitted with a Polaris® valve can undergo an MRI examination, even immediately after the device has been implanted.

The results of in vitro tests have demonstrated that the Polaris® valve does not present any danger in the following examination conditions:

- MRI with a static magnetic field limited to 3-Tesla
- whole body averaged SAR (Specific Absorption Rate) limited to 2.8 W/kg for 15 minutes exposure at 3-Tesla.

The tests conducted in accordance with the standard, ASTM F 2182-02a, showed that the rise in temperature caused by exposure to 3-Tesla MRI was negligible and had no physiological consequences for the patient.

Tests conducted in accordance with standards, ASTM F 2213-06 and ASTM F 2052-06, showed that the torque and displacement force induced by a magnetic field of 3 Tesla or less did not present any risks for the patient.

NOTE:

THE PATIENT MUST BE INFORMED THAT HE/SHE IS LIKELY TO FEEL SLIGHT DISCOMFORT, WHICH IS COMPLETELY INOFFENSIVE, DURING AN MRI EXAMINATION.

Tests have demonstrated that the performance (operating pressures calibration), magnetic properties and functionality (valve can still be adjusted and self-locking mechanism of the rotor remains functional) of the Polaris® valve are not affected by even repeated exposure to MRI examinations MRI from 1.5 Tesla to 3 Tesla (30,000 Gauss).

The Polaris® valve mechanism is designed to prevent accidental changes in operating pressure in the standard conditions for MRI examination at 3 Tesla or less, as long as there is no rotational movement on the valve during MRI exposure.

WARNING:

EVEN THOUGH THE SELF-LOCKING SYSTEM OF THE VALVE WAS DESIGNED TO BE INSENSITIVE TO MAGNETIC FIELDS, IT IS RECOMMENDED THAT THE PRESSURE SETTING OF THE VALVE IS CHECKED BEFORE AND AFTER THE MRI EXAMINATION.

PRECAUTIONS:

DURING MRI EXAMINATIONS, MAKE SURE THE PATIENT REMAINS IMMOBILE IN THE IMMEDIATE PROXIMITY TO THE TUNNEL AND INSIDE IT.

IF THERE IS A ROTATIONAL MOVEMENT ON THE VALVE SIMULTANEOUS WITH EXPOSURE TO A POWERFUL MAGNETIC FIELD (FOR EXAMPLE: 3-TESLA MRI) THERE IS A POSSIBLE RISK OF THE PRESSURE SETTING BEING ACCIDENTALLY CHANGED.

CHOOSE A POLARIS® VALVE IMPLANTATION SITE AWAY FROM AREAS OF SIGNIFICANT CLINICAL INTEREST, SUCH AS A TUMOR, THAT MAY REQUIRE REPEATED FUTURE MRI EXAMINATIONS. INDEED THE MICRO-MAGNETS IN POLARIS® VALVE ARE A POTENTIAL SOURCE OF ARTIFACTS ON MRI IMAGES. THE SIZE OF THESE ARTIFACTS COULD BE VERY LARGE IN SIZE IN RELATION TO THE SIZE AND SHAPE OF THE VALVE.

8. Sterilization – Decontamination of Polaris® Valves and Valve Kits

The Polaris® valves and valve kits are packed individually in double peel-off, sterile, pyrogen-free packaging, sterilized with ethylene oxide.

WARNINGS:

DO NOT USE THE VALVES OR VALVE KITS IF THE STERILE PACKAGING IS OPEN OR DAMAGED, OR IF THE EXPIRATION DATE HAS PASSED.

VALVES AND VALVE KITS ARE SINGLE USE DEVICES. DO NOT RE-STERILIZE OR RE-USE AFTER UNPACKING AND/OR EXPLANATION.

NOTE:

SOPHYSA CANNOT BE HELD RESPONSIBLE FOR THE PERFORMANCE OF ANY PRODUCT THAT HAS BEEN RE-STERILIZED AND/OR RE-USED, NOR FOR ANY COMPLICATIONS WHICH MIGHT RESULT FROM THIS.

9. Cleaning – Decontamination of the Polaris® Adjustment Kit

The various components of the Polaris® adjustment kit are not sterile and must not be sterilized.

However, if there is a low risk of infection, decontamination of the Polaris® adjustment kit may be performed by manual disinfection followed by a final rinse.

Disinfection is best performed by applying a disinfectant, using material (swab) impregnated with a standardized solution of a broad spectrum disinfectant (bactericide, tuberculicide (mycobactericide), fungicide and virucide) in accordance with the recommendations of the manufacturer.

This disinfection must be followed by a meticulous final rinse, again using material (swab), with an appropriate type of water, followed by drying.

The cleaning/disinfection should be performed in accordance with standard hospital procedures.

PRECAUTIONS:

DO NOT STERILIZE OR IMMERSE THE ADJUSTMENT KIT COMPONENTS.

IRREVERSIBLE ALTERATIONS TO THE MARKINGS, DISTORTION OF THE PLASTIC PARTS AND/OR DEMAGNETIZATION OF THE MAGNET RISK RENDERING THE ADJUSTMENT DEVICE UNUSABLE.

ONLY USE STANDARD DISINFECTANT SOLUTIONS IN ACCORDANCE WITH THE METHODS VALIDATED FOR THE DECONTAMINATION OF COMPONENTS OF THE POLARIS® ADJUSTMENT KIT.

NON-ADHERENCE TO STANDARD DECONTAMINATION RULES COULD INDUCE A RISK OF MICROBIAL CONTAMINATION.

10. Instructions

Choice of valve model

The choice of Polaris® adjustable pressure valve model is left to the neurosurgeon and depends upon the clinical needs of the patient.

PRECAUTIONS:

DO NOT USE PRE-CONNECTED VALVE KITS (VALVES WITH PRE-ATTACHED DISTAL CATHETER) FOR VENTRICULO-ATRIAL SHUNTS.

THERE COULD BE IMPLANTATION DIFFICULTIES RELATING TO THE SHUNT LENGTH ADAPTATION AT THE ATRIUM.

USE AN SPVA (ANTECHAMBER) TYPE MODEL, AN SPVB (BURR HOLE RESERVOIR) TYPE MODEL, OR ONE OF THE SPV MODELS COMBINED WITH A VENTRICULAR CATHETER WITH RESERVOIR IF THE PLAN IS TO USE THE SHUNT SYSTEM TO CHECK THE PATENCY OF THE SHUNT, TO SAMPLE THE CSF OR FOR INJECTIONS.

Adjustment of a Polaris® valve BEFORE IMPLANTATION:

selection of operating pressure

Before opening the packaging that guarantees the sterility of the valve, it is necessary to adjust the valve rotor in the position which corresponds to the pressure chosen for the implantation by the surgeon depending upon the needs of the patient.

NOTE:

If a valve is pre-connected to a SiphonX® GRAVITATIONAL ANTI-SIPHON DEVICE,, WHEN ADJUSTING TAKE INTO ACCOUNT THE FACT THAT THE SIPHONX® ADDS UP TO 200mmH₂O TO THE VALVE OPERATING PRESSURE.

The sterile packaging is designed to allow adjustment of the valve through the blister packaging, using a Polaris® adjustment kit.

PRECAUTIONS:

CHANGING THE VALVE PRESSURE MUST ONLY BE DONE BY A NEUROSURGEON, OR OTHER QUALIFIED PERSON.

DO NOT UNPACK OR IMPLANT THE VALVE WITHOUT PREVIOUSLY ADJUSTING IT TO THE OPERATING PRESSURE ESTIMATED TO BE THE MOST APPROPRIATE FOR THE PATIENT.

THE OPERATION SCAR AND/OR POST-OPERATIVE EDEMA CAN MAKE POSITIONING THE LOCATOR PAINFUL AND INACCURATE. THIS COULD THEN CAUSE DIFFICULTIES IN ADJUSTMENT OR EVEN MAKE IT IMPOSSIBLE.

WARNING:

THE ADJUSTMENT KIT MUST NOT BE USED IN THE OPERATIVE FIELD. AS THE ADJUSTMENT KIT CANNOT BE STERILIZED, USING IT DURING THE OPERATION WOULD CAUSE A HIGH RISK OF INFECTION FOR THE PATIENT.

1) Identification of the valve model

Identify the reference and/or the pressure range of the Polaris® valve to be adjusted: Refer to the labeling on the package.


Figure 8: Details of the blister pack of a Polaris® SPV valve set to the highest pressure (for this example)

2) Choice of Locator Reading Area

Use the Locator on which the pressure values for valve models with standard calibration are engraved: 30-200mmH₂O.

For a "special pressure" valve model (10-140, 50-300 or 80-400mmH₂O), fix the reading ring, which corresponds to the valve model, onto the Locator, positioning the two holes in the ring on the two pins on the Locator.

WARNING:


DO NOT USE A POLARIS® ADJUSTMENT KIT WITHOUT PREVIOUSLY IDENTIFYING THE VALVE MODEL AND MAKING SURE THAT THE PRESSURE RANGE SHOWN ON THE LOCATOR CORRESPONDS TO THIS MODEL.

OVER-DRAINAGE OR UNDER-DRAINAGE CAN RESULT FROM TAKING A READING AND/OR MAKING AN ADJUSTMENT WITH A DIFFERENT PRESSURE RANGE ON THE LOCATOR FROM THAT OF THE IMPLANTED VALVE MODEL.

3) Positioning the Locator

Place the **Locator** in the slot provided for it on outside of the packaging tray, centered above the valve, with the green arrow on the Locator pointing in the direction of the CSF flow as shown by the arrow on the valve.


Figure 9

4) Positioning the Compass and Reading the Pressure

Place the **Compass** in the housing provided on the **Locator** using the guide pins: align the lines on the Compass with the pressure values engraved on the Locator. On the **Locator** read the operating pressure shown by the needle on the **Compass**.

Figure 10

PRECAUTION:

WHEN READING THE POSITION OF THE PRESSURE SETTING, MAKE SURE THAT THE MAGNET OR ANY OTHER FERROMAGNETIC OBJECT IS NOT WITHIN A RADIUS OF 0.5M AROUND THE COMPASS SO THAT A FALSE READING IS NOT INDUCED BY THE INFLUENCE OF AMBIENT MAGNETIC FIELDS.

Reading this initial position of the pressure setting can also be performed visually, without using the Polaris® adjustment kit, as the valve rotor is visible through the valve and the blister pack.

5) Positioning the Magnet and Adjusting the Pressure


Remove the **Compass** and insert the **Magnet** into the **Locator**.

Move the green marker on the Magnet to face the initial position of the pressure setting (turn the Magnet in the Locator if necessary).

Hold the **Locator** with one hand.


Slide the **Magnet**, back and forth quickly along the axis of the initial position of the pressure setting, to unlock the valve rotor, then replace it correctly in the center of the **Locator**.

Figure 11

Turn the **Magnet** to the position chosen as the operating pressure for the valve which is to be implanted.

Figure 12

Remove the **Magnet** vertically from the **Locator** to ensure effective re-locking of the valve rotor in the chosen operating pressure position.

NOTE:

DO NOT ATTEMPT TO TURN THE MAGNET IF IT IS NOT CENTERED ON THE LOCATOR, OR IF IT IS NOT COMPLETELY INSERTED INTO THE LOCATOR. UNLOCKING OF THE VALVE ROTOR COULD BE COMPROMISED.

PRECAUTION:


ALWAYS REMOVE THE MAGNET FROM THE LOCATOR VERTICALLY. THIS ENSURES OPTIMUM LOCKING OF THE ROTOR IN ITS NEW POSITION.

FAILURE TO REMOVE THE MAGNET FROM THE LOCATOR VERTICALLY RISKS CHANGING THE SELECTED POSITION AND/OR CAUSING INCORRECT LOCKING OF THE ROTOR.

6) Confirming the operating pressure setting

Checking is performed using the **Compass** in accordance with the procedure described in Step 4.

Figure 13


PRECAUTION:

WHEN READING THE POSITION OF THE PRESSURE SETTING, MAKE SURE THAT THE MAGNET OR ANY OTHER FERROMAGNETIC OBJECT IS NOT WITHIN A RADIUS OF 0.5m AROUND THE COMPASS SO THAT A FALSE READING IS NOT INDUCED BY THE INFLUENCE OF AMBIENT MAGNETIC FIELDS.

Checking can also be performed visually, without using the Polaris® adjustment kit, as the rotor is visible through the valve and the blister pack.

Figure 14

7) Recording the operating pressure chosen for the implantation

Note the pressure value read in Step 6 on the Polaris® Patient Identification Card (PIC).

Implantation Technique

Implantation of a Polaris® valve must take account of current aseptic neurosurgical practices.

The implantation of a shunt including a Polaris® valve may be performed in several ways.

The surgeon will choose the technique depending upon his experience and the clinical status of the patient.

The final implantation of the device must satisfy the conditions for optimal drainage of the CSF.

The surgeon must select the implantation area taking into account the fact that the valve is a potential source of artifacts when an MRI examination is performed (cf. §7 - Behavior during Magnetic Resonance Imaging (MRI)).

PRECAUTION:

DO NOT PERFORM THE IMPLANTATION OF A SHUNT WITHOUT HAVING A REPLACEMENT SHUNT SYSTEM AVAILABLE IN CASE IT IS REQUIRED.

WARNING:

DO NOT CARRY OUT AN ADDITIONAL TEST BEFORE IMPLANTATION: EACH VALVE HAS BEEN INDIVIDUALLY CALIBRATED AND CHECKED. ANY PRE-OPERATIVE PRESSURE TESTS WILL INCREASE THE RISK OF INFECTION.

Ventricular Catheter

- Introduce the catheter into the ventricle using the introducing stylet supplied for this purpose.
- If necessary, adjust the implantation depth of the ventricular catheter with the right angle adapter supplied. Position it in the axis of the burr hole.
- Purge the catheter of air with the CSF.
- If necessary, check that the reservoir is properly filled, and then clamp.
- Connect and delicately ligate the catheter to the inlet connector of the valve (or that of the reservoir for valve models with integrated reservoirs). Check that the arrow located on the upper surface of the valve is correctly oriented in the direction of the flow. The clamp can then be released.

WARNING:

ENSURE THAT THE ARROW ON THE UPPER SURFACE OF THE VALVE IS CORRECTLY ORIENTED IN THE DIRECTION OF THE FLOW: ASSEMBLY OF THE VALVE IN THE OPPOSITE DIRECTION WOULD PREVENT ANY DRAINAGE.

PRECAUTION:

DUE TO THE FRAGILITY OF THE SILICONE, USING METALLIC FORCEPS FOR INSERTING CATHETERS AND FOR LIGATING THEM ONTO CONNECTORS IS NOT RECOMMENDED. THIS WOULD CREATE A RISK OF CUTTING OR PIERCING THE CATHETERS.

Valve

PRECAUTION:

POLARIS® MODELS MUST BE IMPLANTED ON A BONY SURFACE, UNDER SUB-CUTANEOUS TISSUE LESS THAN 8mm THICK.

THE SELECTED IMPLANTATION SITE MUST BE SUFFICIENTLY UNENCUMBERED TO FACILITATE DETECTION OF THE VALVE BY THE LOCATOR WHEN ADJUSTMENTS ARE MADE LATER (SPACE OF AT LEAST 4cm ALL AROUND THE VALVE, AWAY FROM THE EAR AND THE NECK).

IF THESE PRECAUTIONS ARE NOT OBSERVED, THERE MAY BE DIFFICULTIES WITH ADJUSTMENT OR IT MAY EVEN BE IMPOSSIBLE TO MAKE AN ADJUSTMENT, DUE TO POOR LOCALIZATION OF THE VALVE BY THE LOCATOR OR THE MAGNET MOVING THE ROTOR INSUFFICIENTLY.

- Purge the valve of air. To prevent any risk of introducing an air bubble, it is recommended that the valve be left to fill directly with the patient's CSF. In the majority of cases, the valve fills immediately.

However, in patients with low intracranial pressure or if the valve is set to a high pressure, the valve may not fill spontaneously.

In this case:

place a piece of catheter on the outlet connector and slowly aspirate the CSF using a syringe fitted with a Luer connector

or even:

press the dome of the reservoir to allow the CSF to fill the valve (case of SPVA and SPVB type models).

PRECAUTION:

DO NOT FILL OR PURGE THE VALVE WITH ANY LIQUID OTHER THAN THE PATIENT'S CSF OR WATER FOR INJECTION (WFI) BEFORE IMPLANTATION TO AVOID ANY RISK OF DEPOSITS IN THE VALVE, WHICH COULD LEAD TO AN OBSTRUCTION IN THE SHUNT SYSTEM OR A BLOCKAGE IN THE VALVE MECHANISM.

- Check that the valve is correctly filled with CSF and there are no air bubbles inside the valve. If this is not the case, continue to purge. The presence of air bubbles could cause a significant change to the operating pressure initially chosen.
- Check that the arrow on the upper surface of the valve is visible and correctly oriented in the direction of the CSF flow.

PRECAUTION:

DO NOT IMPLANT THE VALVE WITHOUT SUTURING IT TO THE UNDERLYING TISSUES BY ITS TWO CONNECTORS OR BY THE SUTURE HOLES PROVIDED FOR THIS.
IF THE SYSTEM MIGRATES THE DRAINAGE MAY STOP AND OTHER COMPLICATIONS ENSUE.

WARNING:

BEFORE SUTURING THE VALVE TO THE UNDERLYING TISSUES, CHECK THAT THE ARROW ON THE UPPER FACE OF THE VALVE IS VISIBLE.
IF THE ARROW IS NOT VISIBLE THIS MEANS THAT THE VALVE HAS BEEN IMPLANTED THE WRONG WAY UP (UPSIDE — DOWN). IN THIS CASE, IT WILL BE IMPOSSIBLE TO ADJUST THE VALVE AND DIRECT READING WITH THE COMPASS WILL BE REVERSED, FOR EXAMPLE: A HIGH READING (No.5) FOR A VALVE SET TO A LOW POSITION (No.1). THIS COULD CAUSE SERIOUS CLINICAL CONSEQUENCES (OVER- OR UNDER-DRAINAGE).

Case of a valve preconnected to a SiphonX® gravitational anti-siphon device

To prevent any risk of introducing an air bubble, it is recommended that the "Valve + SiphonX®" assembly be left to fill directly with the patient's CSF. For this, make sure it is kept horizontal. In the majority of cases the "valve + anti-siphon device" assembly fills immediately.

PRECAUTION:

FOR OPTIMUM OPERATION OF A POLARIS® VALVE WITH A PRE-ATTACHED SIPHONX® GRAVITATIONAL ANTI-SIPHON DEVICE, MAKE SURE THE ASSEMBLY IS POSITIONED PARALLEL TO THE AXIS OF THE BODY OF THE PATIENT.

Peritoneal Catheter

- Make a short peritoneal incision in the peri-umbilical region.
- Tunnel the distal catheter.
- Connect the proximal end of the catheter to the valve outlet connector and ligate it delicately.
- Check the flow of CSF.
- Adapt the length of the catheter.
- Bury the distal end of the catheter in the peritoneal cavity.

PRECAUTION:

DUE TO THE FRAGILITY OF THE SILICONE, USING METALLIC FORCEPS FOR INSERTING CATHETERS AND FOR LIGATING THEM ONTO CONNECTORS IS NOT RECOMMENDED. THIS WOULD CREATE A RISK OF CUTTING OR PIERCING THE CATHETERS.

Post-operative X-ray examination:

Identification of the Valve Model and Pressure Reading

NOTE:

DURING THE RADIOLOGICAL EXAMINATION ORIENTATE THE PATIENT SO THAT THE RAY SOURCE POINTS PERPENDICULARLY ONTO THE VALVE BODY. IN THIS WAY IDENTIFYING THE VALVE BY ITS RADIO-OPAQUE POINT IS MADE EASY.


The radio-opaque points to the left of the inlet connector are used to identify the pressure range of the valve:

Number of radio-opaque points	zero ○○○	one ○○●	two ○○●●	three ●●●●
Maximum operating pressure (mmH ₂ O)	140	200	300	400

NOTE:


THE PRESSURE RANGE OF THE IMPLANTED VALVE IS SHOWN ON THE PATIENT CARD, IN THE DESCRIPTION OF THE VALVE MODEL USED.

Figure 15: Radiographic identification of the different models and the correlation with the pressure range on the Locator.

It is easy to read the pressures on a radiological image (Figure 16).

Figure 16: Confirmation of the operating pressure (at the top: views of the Polaris® adjustment kit; at the bottom, radiological views of the valve (SPV model))

On the Polaris® valve, each of the five valve pressures is located by the position of the radio-opaque identification point for the corresponding operating pressure.

Locate the valve inlet connector, wider due to the presence of a nut.

The position for the lowest operating pressure (Position No. 1) is the position located nearest to this connector.

If the right-hand valve rotor micro-magnet is facing the radio-opaque point nearest the inlet connector, the valve is set at the lowest pressure (Position No. 1).

Then, moving clockwise away from the inlet connector, each of the following radio-opaque points corresponds to a higher pressure.

So the radio-opaque point furthest away from the valve inlet connector corresponds to the highest pressure (Position No. 5).

Patency Test (post-operative)

There are two steps for the post-operative test on the patency of the shunt:

1. Testing the patency of the ventricular catheter

NOTE:

THIS TEST IS POSSIBLE WITH SPVA (ANTECHAMBER) AND SPVB (BURR HOLE RESERVOIR) TYPE MODELS. FOR THE SPV TYPE MODEL, A VENTRICULAR CATHETER WITH RESERVOIR MUST BE USED.

Pinch the catheter with a finger just after the valve outlet connector.

With another finger, press the reservoir to make the CSF flow back into the ventricular catheter. A reservoir that cannot be compressed easily or does not fill quickly may indicate there is an obstruction in the ventricular catheter.

2. Patency test downstream of the reservoir (valve and distal catheter)

NOTE:

THIS CONTROL IS IMPOSSIBLE WITH A SPVB TYPE MODEL (BURR HOLE RESERVOIR) BECAUSE THERE IS NO ACCESS TO THE VENTRICULAR CATHETER UPSTREAM OF THE RESERVOIR.

Pinch the catheter with a finger just before the reservoir, then with another finger press the reservoir to push the CSF through the valve and distal catheter. A reservoir that cannot be compressed easily may indicate an obstruction either of the valve or the distal catheter.

PRECAUTION:

DO NOT RELY ONLY ON THE CHARACTERISTICS OF THE PATENCY TEST TO DIAGNOSE AN OBSTRUCTION IN THE SHUNT SYSTEM. OBSTRUCTION OF A SHUNT SYSTEM CAN OCCUR IN ANY OF ITS COMPONENTS AND SHOULD BE DIAGNOSED FIRST OF ALL BY THE CLINICAL DATA AND ADDITIONAL EXAMINATIONS.

Sampling the CSF and injection

Access to the CSF is obtained by pricking the reservoir with a 24G (or smaller diameter) Huber needle. The integral reservoir on the SPVA and SPVB type models is designed for occasional use. Its watertight performance is reduced after very frequent pricking into the dome.

- To inject in the proximal direction, compress the catheter just after the valve outlet connector.
- To inject in the distal direction, compress the catheter upstream of the reservoir.

NOTES:

ELECTIVE INJECTION IN THE DISTAL DIRECTION IS NOT POSSIBLE WITH A MODEL OF THE SPVB TYPE (BURR HOLE RESERVOIR) BECAUSE THERE IS NO ACCESS TO THE VENTRICULAR CATHETER UPSTREAM OF THE RESERVOIR.

ENSURE THAT THE BASE OF THE RESERVOIR IS NOT CROSSED WITH THE NEEDLE.

PRECAUTIONS:

DO NOT INJECT INTO, OR TAKE SAMPLES FROM, THE CSF WITHOUT HAVING TESTED THE SHUNT PATENCY.

SIGNIFICANT OVERPRESSURE COULD DAMAGE THE SHUNT IF IT IS KNOWN THAT THERE IS AN OBSTRUCTION.

DO NOT INJECT TOO FAST OR INJECT TOO GREAT A VOLUME.
THE INCREASE IN PRESSURE COULD DAMAGE THE SHUNT.

DO NOT USE A SYRINGE WITH A VOLUME OF LESS THAN 10CC FOR INJECTIONS OR TAKING SAMPLES.
TOO GREAT A PRESSURE COULD DAMAGE THE SHUNT.

Adjustment of a Polaris® valve AFTER IMPLANTATION: reading and/or changing the operating pressure setting

The Polaris® adjustment kit is designed specifically for setting Polaris® valve pressure.

PRECAUTION:

CHANGING THE VALVE PRESSURE MUST ONLY BE DONE BY A NEUROSURGEON, OR OTHER QUALIFIED PERSON.

NOTE:

IF A VALVE IS PRE-CONNECTED TO A SiphonX® GRAVITATIONAL ANTI-SIPHON DEVICE, WHEN ADJUSTING TAKE INTO ACCOUNT THE FACT THAT THE SiphonX® ADDS UP TO 200mmH₂O TO THE VALVE OPERATING PRESSURE.

For the adjustment, the patient should be positioned such that the valve implantation site is easy to access. Having valve horizontal is recommended.

Figure 17

1) Identification of the valve model

Identify the reference and/or the pressure range of the Polaris® valve to be adjusted: Refer to the Patient Identification Card or the patient record (traceability label) and/or perform an X-ray examination of the valve (cf. §10 - Post-operative X-ray examination).

2) Choice of Locator reading area

Use the Locator, by default adapted to valves with standard calibration (30-200mmH₂O).

For a "special pressure" valve model (10-140, 50-300 or 80-400mmH₂O), fix the reading ring corresponding to the pressure range of the valve model identified in Step 1 on the Locator, by positioning the holes in the ring on the two pins on the Locator.

WARNING:

DO NOT USE A POLARIS® ADJUSTMENT KIT WITHOUT PREVIOUSLY IDENTIFYING THE VALVE MODEL AND MAKING SURE THAT THE PRESSURE RANGE SHOWN ON THE LOCATOR CORRESPONDS TO THIS MODEL. OVER-DRAINAGE OR UNDER-DRAINAGE CAN RESULT FROM TAKING A READING AND/OR MAKING AN ADJUSTMENT WITH A DIFFERENT PRESSURE RANGE ON THE LOCATOR FROM THAT OF THE IMPLANTED VALVE MODEL.

3) Positioning the Locator

Palpate the valve implantation site in order to determine both the location and orientation of the valve.


The inlet and outlet connectors at each end of the valve and the reservoir (if fitted) are the easiest components to locate. Locate these first.

Place the **Locator** on the implant site with its axis aligned with that of the valve connectors, and with the green arrow pointing in the direction of the CSF flow.

Center the Locator at best above the valve: Locate the valve by palpating it through the cut-out in the center of the Locator.

Press the Locator firmly against the valve so as to center and immobilize the valve.

Figure 18

4) Positioning the Compass and Reading the Pressure

Hold the Locator with one hand.

Place the **Compass** in the housing provided on the **Locator** using the guide pins: align the lines on the Compass with the pressure values engraved on the Locator. On the **Locator** read the operating pressure shown by the needle on the **Compass**.

Figure 19

PRECAUTION:

WHEN READING THE POSITION OF THE PRESSURE SETTING, MAKE SURE THAT THE MAGNET OR ANY OTHER FERROMAGNETIC OBJECT IS NOT WITHIN A RADIUS OF 0.5M AROUND THE COMPASS SO THAT A FALSE READING IS NOT INDUCED BY THE INFLUENCE OF AMBIENT MAGNETIC FIELDS.

5) Orientation of the Locator

It is recommended that the positioning of Locator in relation to the axis of the valve is fine tuned as follows:

Remove the **Compass** and insert the **Magnet** into the **Locator**.

Move the green marker on the Magnet to face the initial position of the pressure setting (turn the Magnet in the Locator if necessary).

Figure 20

Slide **Magnet** back and forth quickly along the axis of the valve rotor current position as just determined in order to unlock the rotor, then return it to a position centered in the **Locator**.

Figure 21

Turn the **Magnet** slowly, until it goes slightly past the lowest or highest position (No. 1 or No. 5) whichever is furthest away from the initial position.

Figure 22

Remove the Magnet vertically from the Locator, then insert the **Compass** into the **Locator** as described in Step 4.

If the adjustment has been performed correctly, the **Compass** will confirm that the lowest or highest position selected has been reached.

If the Compass needle does not point exactly to this operating pressure value (No. 1 or No. 5), slightly adjust the orientation of the **Locator** (rotation) in order to center the **Compass** needle on the corresponding mark.

If the needle has not moved at all, reattempt the adjustment process after repeating the back and forth motion with the Magnet.

This sequence also makes it possible to confirm that the valve can be adjusted (unlocking and rotation of the rotor).

6) Setting the pressure

Hold the **Locator** firmly in one hand in the precise position defined in Steps 3 and 5, until the last step in the adjustment process is completed.

Remove the **Compass** and insert the **Magnet** into the **Locator**.

Make sure the green marker on the Magnet is aligned with the lowest or highest pressure in Step 5 (turn the Magnet in the Locator if necessary).

Move the **Magnet** quickly backwards and forwards along the axis of the lowest or highest position in Step 5. Then replace it correctly in the center of the **Locator**.

Figure 23

Turn the **Magnet** slowly in the opposite direction until the new position chosen as the operating pressure for the implanted valve is reached.

Figure 24

Remove the **Magnet** vertically from the Locator to ensure effective re-locking of the valve rotor at the new operating pressure position.

NOTE:

DO NOT ATTEMPT TO TURN THE MAGNET IF IT IS NOT CENTERED ON THE LOCATOR, OR IF IT IS NOT COMPLETELY INSERTED INTO THE LOCATOR. UNLOCKING OF THE VALVE ROTOR COULD BE COMPROMISED.

PRECAUTION:

ALWAYS REMOVE THE MAGNET FROM THE LOCATOR VERTICALLY. THIS ENSURES OPTIMUM LOCKING OF THE ROTOR IN ITS NEW POSITION.

FAILURE TO REMOVE THE MAGNET FROM THE LOCATOR VERTICALLY RISKS CHANGING THE SELECTED POSITION AND/OR CAUSING INCORRECT LOCKING OF THE ROTOR.

7) Confirming the operating pressure setting

Checking is performed using the Compass in accordance with the procedure described in Step 4.

Figure 25

Confirmation can also be performed by X-ray.

In the recommended implantation conditions, because direct pressure readings are possible with the adjustment kit, X-ray examination is optional. However, if the information is to be placed in the patient record, an X-ray should be performed.

NOTE:

SETTING THE CHOSEN OPERATING PRESSURE CANNOT ALWAYS BE PERFORMED AT THE FIRST ATTEMPT. IT IS POSSIBLE THAT THE PROCEDURE WILL HAVE TO BE REPEATED ONE OR MORE TIMES, IN ORDER TO ENSURE THAT THE DESIRED PRESSURE SETTING IS OBTAINED.

PRECAUTIONS:

WHEN READING THE POSITION OF THE PRESSURE SETTING, MAKE SURE THAT THE MAGNET OR ANY OTHER FERROMAGNETIC OBJECT IS NOT WITHIN A RADIUS OF 0.5m AROUND THE COMPASS SO THAT A FALSE READING IS NOT INDUCED BY THE INFLUENCE OF AMBIENT MAGNETIC FIELDS.

ADJUSTING THE PRESSURE USING A POLARIS® ADJUSTMENT KIT REQUIRES EXPERIENCE. IF THE USER IS NOT FAMILIAR WITH THE USE OF THE KIT, IT IS RECOMMENDED THAT THE PRESSURE ADJUSTMENT BE CONFIRMED WITH AN X-RAY.

IT IS ADVISABLE TO MONITOR THE PATIENT CAREFULLY FOR 24 HOURS FOLLOWING ANY CHANGE IN THE VALVE OPERATING PRESSURE SETTING.

8) Recording the new operating pressure

Note the pressure value read in Step 7 on the Polarix® Patient Identification Card (PIC).

Potential Causes of Adjustment Difficulties and Troubleshooting

Poor positioning of the Locator in relation to the center of the valve or incorrect orientation of this Locator in relation to the axis of the connectors may lead to an inaccurate or incorrect reading of the operating pressure on the Compass and/or difficulties in making adjustments with the Magnet.

Before anything else, check that the **Locator** is correctly positioned and oriented above the valve (cf. Steps 3 "Positioning the Locator" and 5 "Orientation of the Locator" in the §10 - Adjustment of a Polarix® Valve after implantation).

PRECAUTION:

IT IS NECESSARY TO ENSURE THAT THE LOCATOR IS CORRECTLY CENTERED ABOVE THE VALVE AND THAT ITS BASE IS ABSOLUTELY PARALLEL TO THE VALVE SURFACE (FIGURES 26a & 26b).

BEING OFF-CENTER BY MORE THAN 2MM OR ANY INCLINATION IN RELATION TO THE VALVE COULD CAUSE AN INCORRECT PRESSURE READING OR INCOMPLETE UNLOCKING OF THE ROTOR, MAKING IT IMPOSSIBLE FOR THE ROTOR TO ROTATE AND THUS THE PRESSURE TO BE CHANGED.

Figure 26a



Figure 26b

However, there may be other causes for difficulties in adjustment:

- Excessively thick sub-cutaneous tissue above the valve (>8mm) caused, for example, by post-operative edema, a hematoma or cicatricial tissue would make the localization of the valve inaccurate, thus causing complications in the pressure adjustment procedure.
- The presence of debris or deposits inside the valve could also make adjustment using the usual procedure difficult.

A special procedure may then be applied to facilitate the adjustment of the valve, possibly using a fluoroscopic examination:

- A. Palpate the valve implantation site in order to determine both the location and orientation of the valve.
The inlet and outlet connectors at each end of the valve and the reservoir (if fitted) are the easiest components to locate. Locate these first.
- B. Place the **Compass**, without the Locator, directly onto the implant site, adhering to the orientation of the valve, and centered as well as possible over the valve.
Note the direction of the Compass needle, aligned on the axis of the rotor of the implanted valve.
- C. Remove the **Compass** and place the **Magnet** directly onto the implant site, with the Magnet aligned with the direction of the rotor located in the previous step.

- D. Slide the **Magnet** back and forth quickly along the axis of the current position to unlock the valve rotor. Ensure that the Magnet remains in the plane of the valve. Repeat this procedure several times if necessary until the valve rotor is unlocked. Then replace the **Magnet** in a position centered over the valve.
- E. Then turn the **Magnet** until the desired new operating pressure position is reached. Confirm the new setting using the **Compass**, or with an X-ray.

The implantation of a valve not performed in the recommended conditions (cf. §10 - *Implantation Technique*) may lead to a Compass reading of a pressure value inconsistent with the patient record or the clinical status of the patient.

In this case, an X-ray removes any doubt. It is the absolute proof of a successful adjustment and the correct direction of the implantation.

Finally, pressure adjustment is impossible if the valve is positioned the wrong way up (upside down).

11. Precautions for the Daily Life of the Patient

A Patient Identification Card is supplied with the Polaris® valve. It enables the neurosurgeon to consult and update information relating to the implanted device (reference, operating pressure, implantation site, etc.) systematically and to ensure that the illness is properly monitored.

PRECAUTION:

THE PATIENT SHOULD BE WARNED THAT IT IS IMPORTANT TO CARRY HIS/HER PATIENT IDENTIFICATION CARD (PIC) AT ALL TIMES.

THIS CARD GIVES INFORMATION ON THE MEDICAL SITUATION OF THE PATIENT TO ALL MEDICAL PERSONNEL.

The magnetic self-locking system of the Polaris® valve is designed to make the magnetic rotor insensitive to the influence of standard magnetic fields.

As a result, the following are not likely to affect the valve operating pressure:

- magnetic fields generated by walk-through scanners in airports, microwave ovens, cordless telephones, high tension cables, and TV
- permanent household magnets such as those present in toys, audio headsets and loudspeakers
- magnetic fields created by electric motors operating in equipment such as razors, hairdryers, lawn mowers, etc.

The doctor is responsible for informing the patient or his/her family that the person fitted with a shunt must avoid any activity that may subject this shunt to direct shocks (violent sports, etc.) as these are likely to damage it.

PRECAUTION:

THE PATIENT MUST BE WARNED THAT VIBRATIONS DUE TO THE CSF FLOW MAY POSSIBLY BE FELT BECAUSE OF THE IMPLANTATION OF THE VALVE ON THE SKULL.

12. Complications / Side effects

Complications which may result from the implantation of a CSF shunt system include the inherent risks in the use of drugs, any surgical intervention and the insertion of a foreign body.

PRECAUTION:

PATIENTS TREATED WITH A SHUNT SYSTEM MUST BE CLOSELY MONITORED POST-OPERATIVELY IN ORDER TO DETECT ANY SIGNS OF COMPLICATIONS EARLY.

The doctor is responsible for educating the patient or his/her family about CSF shunt systems, in particular describing the complications linked to implanted shunt systems as well as giving explanations about possible alternative therapies.

The main complications of shunts are obstruction, infection and over-drainage. These complications require rapid intervention by the doctor.

Obstruction

Obstruction is the most frequent complication in shunt systems. It can occur at any point in a shunt.

The ventricular catheter can be obstructed by a blood clot, cerebral tissue or even tumoral cells.

The end of the ventricular catheter can also become embedded in the choroid plexus or in the ventricular wall, either directly or following a collapse of the walls, a consequence of over-drainage.

The cardiac catheter can be colonized by a thrombus while the appearance of a clot around the catheter could cause an embolism in the pulmonary circulation.

The peritoneal catheter may become obstructed by the peritoneum or by intestinal loops.

Loss of patency in a shunt may also be the result of an obstruction by fragments of cerebral tissue or by biological deposits (protein deposits, etc.)

Obstruction of the shunt will quickly result in the reappearance of the signs and symptoms of intracranial hypertension.

These signs and symptoms vary from patient to patient and over time.

In infants and young children, the symptoms may be an abnormal increase in the size of the skull, a bulge in the fontanelles, dilation of the scalp veins, vomiting, irritability with a lack of attention, downward deviation of the eyes, and sometimes convulsions.

In older children and adults, intracranial hypertension due to hydrocephalus may be the cause of headaches, vomiting, blurred vision, diplopia, drowsiness, slowing of movements, gait disorders or psychomotor slowing which could lead to total invalidity.

If an obstruction is confirmed and a patency test does not make it possible to reduce the obstruction, revision surgery or removal of the device must be envisaged.

Infection

Chronic malfunction of the shunt could cause a leak and a discharge of CSF along its length increasing the risk of infection.

Local or systemic infection is another possible complication of CSF shunt systems. It is generally secondary to the colonization of the shunt by cutaneous germs. Nevertheless, as for all foreign bodies, any local or systemic infection can colonize the shunt. Erythema, edema and skin erosions along the length of the shunt may be an indication of an infection of the shunt system.

Prolonged, unexplained fever may also be the result of a shunt system infection.

Septicemia, favored by an alteration in general status, can start from a shunt infection.

If there is infection, removal of the system is indicated in conjunction with the start of a specific treatment by a general or intrathecal route.

Overdrainage

Overdrainage can result in a collapse of the ventricles (slit ventricle syndrome) and the appearance of a subdural hematoma.

In children, depression of the fontanelles, overlapping of the skull bones, even a craniostenosis or a change from communicating hydrocephalus to obstructive hydrocephalus by stenosis of the Aqueduct of Sylvius could occur.

Adults can present with a variety of symptoms such as vomiting, auditory or visual disorders, drowsiness or even headaches in the upright position but which improve in the supine position.

Depending on clinical observations and medical imaging, the doctor can reduce the symptoms of overdrainage and correct the ventricle size by changing the operating pressure of the Polaris® valve. However, immediate drainage of a subdural hematoma may be indicated.

Other

Failure of a shunt system may also be linked to disconnection of its various components.

The ventricular catheter may migrate inside the ventricle. The peritoneal catheter may migrate into the peritoneal cavity under the action of the peristaltic waves of the intestine, while an atrial catheter may migrate into the right-hand cavities of the heart following the blood flow.

Perforation or occlusion of abdominal viscera by the peritoneal catheter could occur.

Growth of the body may progressively cause the catheters to exit their insertion sites.

These malfunctions require the shunt to be replaced immediately.

Cases of cutaneous necrosis over the implantation site are possible.

Over time fibrous adhesions may fix the ventricular catheter in the choroid plexus or the cerebral tissue. If removal is envisaged, gentle rotation of the catheter about its axis may make it possible to free it. The catheter should never be withdrawn forcibly. If it cannot be taken out without forcing, it is preferable to leave it in place rather than risk an intra-ventricular hemorrhage.

Cases of allergy to silicone have been described.

Cases of epilepsy after implantation of a ventricular shunt have been described.

The ruby ball in the valve can potentially take up an off-center position on its housing due to the presence of an aggregation of cells or a protein deposit. Among others, such situations can cause:

- loss of regulatory function in the valve potentially increasing the risk of overdrainage.
- an impaired anti-reflux function.

The mobility of the rotor could be impeded by an aggregation of cells or a protein deposit. This could make it impossible to adjust the valve with the Magnet.

Blood clots, cerebral cells or tumoral cells contained in the CSF could lodge in the valve mechanism, which would have the potential to cause changes in the operating characteristics of the valve.

13. Guarantee

The performance of the Polaris® adjustment kit is only guaranteed with the range of Polaris® valves, designed tested and manufactured by Sophysa.

Sophysa guarantees that this medical device is free from defects in material or manufacturing. Apart from this guarantee, Sophysa does not grant any other guarantee, express or implicit, including marketing or adaptation for a specific use. Sophysa cannot be held responsible for any incident, complication, damage or prejudice resulting directly or indirectly from the use of this device. Sophysa does not authorize anyone to take responsibility on its behalf for its products.

The performance of Polaris® valves is only guaranteed with the range of silicone catheters and accessories designed, tested and manufactured by Sophysa. However, it is possible to use other brands of catheters provided that their internal diameter is identical to that of the catheters recommended by Sophysa.

14. Processing of the products after use

Destruction after use

An unpacked, used or explanted Polaris® valve must be destroyed in accordance with the procedures in force in the medical establishment.


Return of products

If an explanted valve needs to be returned to Sophysa for analysis, it must be returned immersed in water (non-saline), indicating if necessary whether cleaning has been performed.

Never use a saline solution likely to form deposits in the valve body which could block the rotor.

In order to properly assess the returned product, it must be accompanied by an explanatory Return to Manufacturer Authorization form.

15. Symbols

	Catalog reference
	Manufacturer
	Caution, see the Instructions for Use
	Sterilization Method using Ethylene Oxide
	Do not re-use
	Do not re-sterilize
	Use until
	Batch code
SN	Serial number
	CE Conformity Marking
	Powerful magnet

References

Name/Description	Reference
POLARIS® ADJUSTABLE PRESSURE VALVES FOR CSF SHUNTING	
Polaris® adjustable valve; 30-200 mmH₂O Adjustable pressure valve, 5 pressures: 30 (Low), 70, 110 (Medium), 150, 200 (High) mmH ₂ O.	SPV
Polaris® adjustable valve; 10-140 mmH₂O Adjustable pressure valve, 5 pressures: 10 (Low), 40, 80 (Medium), 110, 140 (High) mmH ₂ O.	SPV-140
Polaris® adjustable valve; 50-300 mmH₂O Adjustable pressure valve, 5 pressures: 50 (Low), 100, 150 (Medium), 220, 300 (High) mmH ₂ O.	SPV-300
Polaris® adjustable valve; 80-400 mmH₂O Adjustable pressure valve, 5 pressures: 80 (Low), 150, 230 (Medium), 330, 400 (High) mmH ₂ O.	SPV-400
Polaris® adjustable valve SPV / Antechamber SPV Valve, 5 pressures: 30 (low), 70, 110 (Medium), 150, 200 (High) mmH ₂ O, with integrated antechamber.	SPVA
Polaris® adjustable valve SPV-140 / Antechamber SPV-140 Valve, 5 pressures: 10 (Low), 40, 80 (Medium), 110, 140 (High) mmH ₂ O, with integrated antechamber.	SPVA-140
Polaris® adjustable valve SPV-300 / Antechamber SPV-300 Valve, 5 pressures: 50 (Low), 100, 150 (Medium), 220, 300 (High) mmH ₂ O, with integrated antechamber.	SPVA-300
Polaris® adjustable valve SPV-400 / Antechamber SPV-400 Valve, 5 pressures: 80 (Low), 150, 230 (Medium), 330, 400 (High) mmH ₂ O, with integrated antechamber.	SPVA-400
Polaris® adjustable valve SPV / Burr hole reservoir SPV Valve, 5 pressures: 30 (Low), 70, 110 (Medium), 150, 200 (High) mmH ₂ O, with integrated burr hole reservoir.	SPVB
POLARIS® ADJUSTABLE PRESSURE VALVES FOR CSF SHUNTING WITH PRE-ATTACHED SIPHONX® GRAVITATIONAL ANTI-SIPHON DEVICE	
Polaris® adjustable valve, 30-200 with SiphonX® SPV Valve, 5 pressures: 30 (Low), 70, 110 (Medium), 150, 200 (High) mmH ₂ O with pre-attached 200mmH ₂ O gravitational anti-siphon device.	SPV-SX
Polaris® adjustable valve, 30-200 with antechamber and SiphonX® SPV valve, 5 pressures: 30 (Low), 70, 110 (Medium), 150, 200 (High) mmH ₂ O, with integrated antechamber and pre-attached 200mmH ₂ O gravitational anti-siphon device.	SPVA-SX
Polaris® adjustable valve, 30-200 with burr hole reservoir and SiphonX® SPV valve, 5 pressures: 30 (Low), 70, 110 (Medium), 150, 200 (High) mmH ₂ O, with integrated burr hole reservoir and pre-attached 200mmH ₂ O gravitational anti-siphon device.	SPVB-SX
Polaris® adjustable valve, 10-140 with SiphonX® SPV Valve, 5 pressures: 10 (Low), 40, 80 (Medium), 110, 140 (High) mmH ₂ O with pre-attached 200mmH ₂ O gravitational anti-siphon device.	SPV-140-SX
Polaris® adjustable valve, 10-140 with antechamber and SiphonX® SPV valve, 5 pressures: 10 (Low), 40, 80 (Medium), 110, 140 (High) mmH ₂ O, with integrated antechamber and pre-attached 200mmH ₂ O gravitational anti-siphon device.	SPVA-140-SX
SHUNT KITS FOR CRANIAL IMPLANTATION	
Complete adjustable Polaris® kit; 30-200mmH₂O SPV Valve, 5 pressures: 30 (Low), 70, 110 (Medium), 150, 200 (High) mmH ₂ O with Pre-attached peritoneal catheter (B905S). Right ventricular catheter (B019-10)	SPV-2010
Complete adjustable Polaris® kit; 30-200mmH₂O / Antechamber SPVA Valve, 5 pressures: 30 (Low), 70, 110 (Medium), 150, 200 (High) mmH ₂ O with Pre-attached peritoneal catheter (B905S). Right ventricular catheter (B019-10)	SPVA-2010
Complete adjustable Polaris® kit; 30-200mmH₂O / Burr hole reservoir SPVB Valve, 5 pressures: 30 (Low), 70, 110 (Medium), 150, 200 (High) mmH ₂ O with Pre-attached peritoneal catheter (B905S). Right ventricular catheter (B019-10)	SPVB-2010
DEVICE FOR READING AND ADJUSTING THE OPERATING PRESSURE OF THE POLARIS® VALVE	
Polaris® Locating Instrument (Locator)	PAK-LI
Polaris® Setting Instrument (Magnet)	PAK-SI
Polaris® Reading Instrument (Compass)	PAK-RI
Polaris® valve adjustment kit (including the three previous components)	PAK

Technical specifications and List of product references may be modified without notice.
Availability may vary according to country.

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Sophysa USA Inc.: