

# SyvekNT™

## Provides Superior Control of Bleeding Wounds for Percutaneous Procedures



### SyvekNT™ is Non-Immunogenic with No Known Contraindications

SyvekNT utilizes poly-N-acetyl glucosamine (pGlcNAc) fibers derived from microalgae, proven to accelerate platelet activation, RBC aggregation and activation, and vasoconstriction.<sup>1</sup>

- Proven results in anticoagulated patients<sup>2</sup>
- Reduced time to ambulation<sup>3</sup>
- No foreign materials remain at puncture site<sup>4</sup>

### Effectively Controls Bleeding for Percutaneous Catheterization

SyvekNT effectively controls bleeding for patients undergoing diagnostic and interventional procedures including femoral and radial catheterization, and is indicated for patients on anticoagulation therapy.

### Applicable for Cardiac Catheterization Procedures

In a published U.S. clinical evaluation of SyvekNT for 1,000 consecutive patients:<sup>2</sup>

- Included stenting, PTCA, electrophysiology and diagnostic catheterization procedures
- Catheter sheaths ranged from 4-12 French
- 76% of interventional patients sheaths were pulled with ACTs > 200
- < 1.5% minor complications

### Can Safely Reduce Time to Ambulation


Clinical Evidence: In a U.S. study of 200 patients undergoing a diagnostic coronary angiography using 6 French catheters:<sup>3</sup>

- 98% were able to safely ambulate after one hour of bedrest
- No major adverse events and only
- 2% minor adverse events\*
- Shortened bedrest may improve patient satisfaction
- May improve catheterization laboratory efficiency

Visit [www.syvek.com](http://www.syvek.com) or call 1.888.666.2560 to learn more about the effectiveness of SyvekNT™.

# Easy to Use and Apply

SyvekNT™ is a soft, white, sterile non-woven pad of poly-N-acetyl glucosamine fibers isolated from microalgae. The 3 cm x 3 cm dressing is attached to medical-grade foam backing and stored in a sterile, single-use package.



**1** Open single pack onto sterile field.



**2** Apply SyvekNT and pressure to the site.

## Instructions for Use

- Open SyvekNT™ and deliver it to the sterile field
- Hold occlusive pressure proximal to the puncture site and remove the sheath
- Allow 1-2 drops of blood from site
- Immediately place SyvekNT over the puncture site and hold firm pressure, while continuing to hold proximal pressure
- Slowly decrease amount of proximal pressure after 4 to 6 minutes (3 – 4 minutes for diagnostic procedures) while continuing to apply firm compression over the puncture site
- Hold firm pressure on the puncture site until bleeding is controlled, then slowly release
- Assess the site without disturbing SyvekNT, then cover with a transparent dressing

## Removing SyvekNT™

- Soak SyvekNT with water
- Gently peel off SyvekNT

\*All of the minor adverse events were managed successfully with additional bedrest of 1 - 2 hours.<sup>3</sup>

<sup>1</sup>Valeri C, Vournakis J. *J Trauma*. 2011;71:S162-S166.

<sup>2</sup>Nader R, Garcia J, Drushal K, et al. *J Invas Cardiol*. 2002; 14:305-307

<sup>3</sup>Palmer B, Gantt D, Lawrence M, et al. *Am J Cardiol*. 2004; 93:96-97

<sup>4</sup>Hirsch JA, Reddy SA, Capasso WE, et al. *Tech Vasc Interv Radiol*. 2003;6:92-95.

## Indications for Use

SyvekNT is intended for the promotion of rapid control of bleeding in patients following hemodialysis and in patients on anticoagulation therapy. SyvekNT is intended for use in the local management of bleeding wounds such as vascular access site, percutaneous catheters or tubes and surgical debridement.

## Safety Information

PRECAUTIONS: There are no known contraindications to the use of SyvekNT. Longer compression may be necessary for hypertensive or obese patients.

## LATEX FREE

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner

© 2012 Marine Polymer Technologies, Inc. All Rights Reserved. Marine Polymer Technologies, Inc., the M logo, mrdh, Syvek, Talymed and MP719 are trademarks or registered trademarks of Marine Polymer Technologies, Inc. in many jurisdictions worldwide. Other product or service names may be trademarks of Marine Polymer Technologies, Inc. or other companies.

## SyvekNT™ Ordering Information

Item #	Quantity	Description
400-09	1 Case of 30 (6 Boxes of 5)	SyvekNT™ 3 cm x 3 cm

Order Fax: 877.270.8500

Order Telephone: 888.666.2560

## About Marine Polymer Technologies, Inc.

SyvekNT™ is manufactured by Marine Polymer Technologies, Inc., an innovative medical device company focused on the development and delivery of biomaterials for the medical field. Since 1992, we have introduced a variety of groundbreaking products in the areas of hemostasis and wound healing to support healthcare providers and their patients. Our evolutionary products are based on the isolation, characterization and formulation of our novel and proprietary polysaccharide polymer (pGlcNAc). We are committed to the highest quality standards that meet both our customers' expectations and regulatory requirements. Marine Polymer Technologies, Inc. complies with FDA QS 820 regulations, the EC Directive 93/42/EEC and is ISO 13485 certified.

## SyvekNT™

Marine Polymer Technologies, Inc.  
107 Water Street  
Danvers, MA 01923  
888.666.2560 | [information@syvek.com](mailto:information@syvek.com)