

## DIRECTIONS FOR APPLICATION AFTER CARDIAC CATHETERIZATION PROCEDURES



1. Open Syvek<sup>®</sup>Patch**SP** pouch and deliver it to the sterile field.
2. Hold occlusive pressure proximal to the puncture site and remove the sheath.
3. Allow 1-2 drops of blood from site.
4. Immediately place Syvek<sup>®</sup>Patch**SP** (tan side AWAY from skin) over the puncture site and hold firm pressure on the Syvek<sup>®</sup>Patch**SP**, while continuing to hold proximal pressure.
5. Slowly release occlusive proximal pressure after 4 to 6 minutes (3 – 4 minutes for diagnostic procedures) while continuing to apply firm compression over the puncture site.
6. Hold firm pressure on the puncture site until hemostasis is achieved, then slowly release.
7. Assess the site without disturbing the Syvek<sup>®</sup>Patch**SP**, then cover with a transparent dressing.

REMOVE Syvek<sup>®</sup>Patch**SP** after 24 Hours:

1. Soak Syvek<sup>®</sup>Patch**SP** with water.
2. Gently peel off Syvek<sup>®</sup>Patch**SP**.

## DESCRIPTION

Syvek<sup>®</sup> Patch**SP** is a soft white, sterile non-woven pad of a cellulosic polymer (poly-N-acetylglucosamine) isolated from a microalgae, attached to medical grade foam backing.

## INDICATIONS

Syvek<sup>®</sup> Patch**SP** is intended for use under the direction of a healthcare professional.

Syvek<sup>®</sup> Patch**SP** promotes the rapid control of bleeding in patients following hemodialysis and in patients on anticoagulation therapy.

Syvek<sup>®</sup> Patch**SP** is indicated for use in the local management of bleeding wounds such as vascular access site, percutaneous catheters or tubes and surgical debridement.

Store at controlled room temperature  
15°C to 30°C (59°F to 86°F).

Sterilized with radiation.

Do not resterilize and/or reuse.  
One time use only.

## PRECAUTIONS

There are no known contraindications to the use of Syvek<sup>®</sup> Patch**SP**.

Longer compression may be necessary for hypertensive or obese patients.

LATEX FREE

## DISTRIBUTED BY



MPS 500-25 rev 04



## DIRECTIONS FOR APPLICATION AFTER CARDIAC CATHETERIZATION PROCEDURES