

MAINTENANCE AND WARRANTY GUIDE

AQUAROLL® A2S1 SYSTEM



INITIAL INSTALLATION

Professional installation of the AquaRoll® A2S1 system is highly recommended. Service and installation contracts are available, upon request.

- **Caution:** Glass beads must be thoroughly rinsed three (3) times with one (1) gallon of water for one (1) minute per rinse prior to initial client use of the AquaRoll® A2S1 system.
- **Caution:** Only use glass beads supplied by AquaRoll®, as they are specially designed for full functionality.
- **Caution:** Do not exceed two and a half (2.5) bags of glass beads at one time in the system.
- **Note:** Fewer beads may be used when treating smaller extremities. The maximum quantity of beads may be utilized when treating large feet, such as with male athletes.

STORAGE REQUIREMENTS

Always press the “On/Off” power switch to “Off” and unplug the device after each daily usage session.

Store the AquaRoll® A2S1 system in a moderate, dry location.

- **Caution:** Never store the AquaRoll® A2S1 system in a location with a temperature below fifty (50) degrees Fahrenheit.
- **Caution:** Never allow unattended children or mentally disabled clients around the AquaRoll® A2S1 system, as there is a risk of ingestion of glass beads.

It is recommended to cover the device, when not in use.

Where necessary to discard the AquaRoll® A2S1 system, ensure full compliance with local waste regulations.

MAINTENANCE

The drainage valve of the tank must be cleaned quarterly, at a minimum. Refer to the AquaRoll® Care and Cleaning Instructions (90000007) for further details.

Glass beads occasionally require replacement after several thousand treatment sessions. Consult an AquaRoll® service professional for further details.

- **Caution:** Failure to replace the glass beads can result in decreased client benefit, due to the loss of uniformity and constancy of the percussive hydro-massage. Such failure can also result in worn-down beads getting stuck inside the drainage and system tubing and can void the warranty.

WARRANTY INFORMATION

The AquaRoll® A2S1 system comes with a one (1) year warranty on parts.

The warranty is void upon the following conditions:

- **improper use and/or evidence thereof**
- **attempted opening and/or repair of the machine by anyone other than the designated service professional**

Additional parts and labor warranties are available for purchase, upon request. Please contact AM Rejuvenation for details.

AM Rejuvenation, LLC

6538 Collins Avenue, #313
Miami Beach, FL 33141

To order or for info, call:
888.850.3797

Find more info online:
www.aquarollUSA.com

Request info by email:
sales@aquarollUSA.com