

Instructions for Use:

Note: Always wear appropriate personal protective clothing as indicated by Good Laboratory Practice. If blood is spilled follow your laboratory policy or national guidance, for safe disposal.

1. Mix the blood sample thoroughly but gently, by at least 8 full inversions of the container.
2. **Pre-citrated blood sample:** using a transfer pipette add 1.25ml of the well mixed sample into the filling reservoir. Ensure the sample volume does not go beyond the filling line (diagram 1).
3. **EDTA blood sample:** In a 12 x 75mm test tube add 4 parts of the blood sample to 1 part Normal saline solution (0.85%-0.145 mol/L, NaCl). Mix very thoroughly, but carefully, by at least 8 full inversions of the test tube. Transfer 1.25ml of this mixture to the blue filling reservoir – ensure that the volume of sample reaches but does not go beyond the filling line (diagram 1).
4. Without delay gently insert the Dispette into the reservoir and observe the blood sample being pushed up inside the bore until it reaches to or beyond the zero level (see diagrams 2 & 3). The Dispette **must** be fully inserted to the bottom of the reservoir – if the blood level in the Dispette then reaches above the zero mark on the graduated scale, apply gentle pressure with the thumb to the side of the tube (see diagram 4) and the blood level will drop – stop at exactly the zero mark.
5. Place the full Dispette assembly (Dispette plus reservoir) in a correctly levelled ESR stand ensuring that the Dispette is at 90 degrees, plus or minus no more than one degree to the horizontal, (i.e. vertical, see diagrams 5a & 5b).
Immediately start the timer.
6. At the end of the timed hour, the result is read by aligning the eye to the level where the red cell column has dropped – leaving clear plasma above (see illustration d on page 3). Record the number of mm the red cells have dropped from the mm scale printed on the Dispette.
7. Record the result as ESR = x mm
or Sed-Rate = x mm

Note: It is common, though an incorrect practice, to record ESR results as x mm/hour or x mm per hour. Both the ICSH and CLSI point out that despite the common title Erythrocyte Sedimentation *Rate*, it is not in fact a rate but a measurement at one hour. Hence both these official organizations currently recommend that results should be reported as: ESR = x mm.




Diagram 1.




Diagram 2.




Diagram 3.






Diagram 4.



**Diagrams
5a & 5b**



PRECAUTIONS:

- a) Thorough mixing of the blood sample before adding to the reservoir is essential. However it must be done gently, shaking the sample will result in haemolysis which may obscure the end point.
- b) Shaking the sample may also create bubbles which will seriously affect the result. If bubbles are present in the Dispette repeat the test with a fresh sample that has been more carefully mixed.
- c) The ESR must be performed at room temperature (defined by the ICSH as 18 to 25°C or 64 to 77°F). Do not place the stand near a window, in direct sunlight or where it may be subjected to drafts.

Note: Should it be impossible to perform the test within this temperature range request a temperature correction chart from Arkray USA Inc.

- d) The ESR is affected by vibration; ensure the stand is placed well away from machinery and that the bench is not subjected to knocks. Remember that vibration may only occur intermittently and/or may be as a result of machinery (e.g. a centrifuge) sited in another part of the building.
- e) Do not pick up the stand to read the result as this will affect other tests in progress. Bring the eye to the level of the red cells to read accurately from the scale (see illustration d).
- f) If the result falls between 130 and 155mm, the exact end point may be obscured by the filling reservoir. In this case report the result as >130 mm. If a numerical value is requested estimations can be made, by reading from the scale of an unused Dispette pipette held next to the test assembly.
Note: Before reporting an estimated result be sure to verify that this complies with your laboratory's best practice and standards. Record and report that the result has been estimated.
- g) Occasionally the level of the red cells is not clear-cut and a 'Christmas Tree' effect may be observed. In such cases the level where the red cells become fully concentrated should be recorded (see illustration e).
- h) In cases of serious infection or leukemia a heavy layer of white cells may be present on top of the column of red cells. This should be ignored and the reading taken from the level of the red cells only.
- i) It must be remembered that the red cells will continue to settle after one hour has elapsed. Therefore it is most important to read the test at exactly one hour after setting up the test in the stand (e.g. the test reading will generally be higher at 1 hour 15 minutes than at 1 hour).
- j) **Note:** The rate of fall of the red cells is not linear. Therefore it must not be assumed that the reading at 60 minutes will be twice that at 30 minutes (see Physical basis of blood sedimentation on page 1).




Illustration (d).

Bring the eye to the level of the red cells to correctly read the result (see paragraph e)

Please note: Illustrations are not to scale

Illustration (e).

Reading the maximum concentration of red cells where 'Christmas Tree' effect is present (see paragraph g).

