
BD Vacutainer®

Rapid Serum Tube (RST)

Blood Collection Tube

REF 368774

Code	Language	Pages	Code	Language	Pages
EN	English	2 - 4	KO	한국어	47 - 49
BG	Български	5 - 7	LT	Lietuvių k.	50 - 52
CS	Česky	8 - 10	LV	Latviski	53 - 55
DA	Dansk	11 - 13	NL	Nederlands	56 - 58
DE	Deutsch	14 - 16	NO	Norsk	59 - 61
EL	Ελληνικά	17 - 19	PL	Polski	62 - 64
ES	Español	20 - 22	PT	Português	65 - 67
ET	Eesti	23 - 25	RO	Română	68 - 70
FI	Suomi	26 - 28	RU	Русский	71 - 73
FR	Français	29 - 31	SK	Slovensky	74 - 76
HR	Hrvatski	32 - 34	SL	Slovenski jezik	77 - 79
HU	Magyar	35 - 37	SV	Svenska	80 - 82
IS	Íslenska	38 - 40	TR	Türkçe	83 - 85
IT	Italiano	41 - 43	UK	Українська	86 - 88
KK	Қазақша	44 - 46	ZH	简体中文	89 - 91
			Symbol and Mark Key		92

INTENDED USE

The BD Vacutainer® Rapid Serum Tube (RST) Blood Collection Tube is a single use tube used to collect, separate, transport and process venous blood specimens to obtain serum for chemistry determinations for *in vitro* diagnostic use. It is used in settings where a venous blood sample is collected by a trained healthcare worker.

Clot time has not been established for patients on heparin therapy, direct thrombin inhibitor therapy or with Factor I deficiency. As a result, the use of the BD Vacutainer® RST Tube in these patients is not recommended.

SUMMARY AND EXPLANATION

The BD Vacutainer® RST Tube is a 13x100mm, 5.0 mL, plastic evacuated tube with an orange BD Hemogard™ Closure. The interior of the tube wall is coated with thrombin to promote rapid clotting. The tube contains a gel barrier polymer at the tube bottom. The density of this material causes it to move upward during centrifugation to the serum-clot interface, where it forms a barrier separating serum from the clot. Serum may be aspirated directly from the collection tube, eliminating the need to transfer to another container. The tube interior is sterile and the product is for single use only. See Instructions For Use, Analytic Equivalency, Limitations of System, Cautions and Warnings.

STORAGE

Store BD Vacutainer® RST Tubes at 4 - 25°C (39 - 77°F). Do not use tubes after their expiration date.

INSTRUCTIONS FOR USE

Open barrier bag using the tear notch on either right or left side of bag.

1. Prepare venipuncture site with an appropriate antiseptic. Use your institution's recommended procedure for standard venipuncture technique and sample collection.
2. Remove needle shield.
3. Perform venipuncture.
4. Place tube in holder and push tube forward until tube stopper has been penetrated.
5. Release tourniquet as soon as blood appears in tube.
6. Wait until tube has filled to its stated volume and blood flow ceases.
7. Pull tube off needle inside holder.
8. Remove tube from holder.
9. Mix **immediately** and **gently** 5 to 6 times by inverting the tube. One inversion is turning the tube upside-down and returning it to its upright position.
10. Allow blood to clot for 5 minutes between 23° and 27°C (73-81°F).
11. Centrifuge tubes at 1500 – 2000g for 10 minutes, at 2000g for 4 minutes, or at 4000g for 3 minutes between 23° and 27°C (73-81°F).

Tubes may also be centrifuged at alternate conditions as validated by the laboratory.

Note: Always follow your facility's protocol for order of draw.

INSTRUCTIONS FOR REMOVAL OF BD HEMOGARD™ CLOSURE

1. Grasp the tube with one hand, placing the thumb under the BD Hemogard™ Closure (for added stability, place arm on a solid surface). With the other hand, twist the BD Hemogard™ Closure while simultaneously pushing up with the thumb of the other hand **ONLY UNTIL THE TUBE STOPPER IS LOOSENED**.
2. Move thumb away before lifting closure. **DO NOT** use thumb to push closure off the tube. **If the tube contains blood, an exposure hazard exists.** To help prevent injury during removal, it is important that the thumb used to push upward on the closure be removed from contact with the tube as soon as the BD Hemogard™ Closure is loosened.
3. Lift closure off tube. In the unlikely event of the plastic shield separating from the rubber stopper, **DO NOT REASSEMBLE CLOSURE**. Carefully remove rubber stopper from tube.

INSTRUCTIONS FOR REINSERTION OF BD HEMOGARD™ CLOSURE

1. Place closure over tube.
2. Twist and push down until stopper is fully resealed. Complete reinsertion of stopper is necessary for the closure to remain securely on the tube during handling.

LIMITATIONS OF SYSTEM

The quantity of blood drawn varies with altitude, ambient temperature, barometric pressure, tube age, venous pressure, and filling technique.

BD Vacutainer® RST Tubes are not recommended for amino acid determinations. Therapeutic drug monitoring (TDM), blood banking and infectious disease performance has not been established, except for anti-CMV IgG and anti-CMV IgM.

Clot time has not been established for patients on heparin therapy, direct thrombin inhibitor therapy or with Factor I deficiency. As a result, the use of the BD Vacutainer® RST Tube in these patients is not recommended.

The flow properties of the barrier material are temperature-related. Centrifuge tubes between 23° and 27°C (73-81°F).

Tubes should not be re-centrifuged once barrier has formed.

CAUTIONS AND WARNINGS

Cautions

1. Do not use tubes if foreign matter is present.
2. Since BD Vacutainer® RST Tubes contain chemical additives, it is important to avoid possible backflow from the tube during blood collection, to minimize the risk of adverse patient reactions. To guard against backflow, observe the following precautions:
 - a. Place patient's arm in a downward position.
 - b. Hold tube with the stopper uppermost.
 - c. Release tourniquet as soon as blood appears in tube.
3. Do not shake. Vigorous mixing may cause foaming or hemolysis.
4. If tubes are not mixed 5 to 6 times immediately after collection, incomplete separation of serum may occur. This may also result in delayed clotting and fibrin formation.
5. Separation of serum from cells by centrifugation should take place within 2 hours of collection to prevent erroneous test results.
6. Remove stoppers with a twist and pull motion. Removal by rolling with the thumb is not recommended.
7. After venipuncture, the top of the stopper may contain residual blood. Take proper precautions when handling tubes to avoid contact with this blood.
8. Overfilling or underfilling of tubes will result in an incorrect blood-to-additive ratio and may lead to incorrect analytic results or poor product performance.
9. Always use appropriate centrifuge carriers or inserts. Use of tubes with cracks or excessive centrifugation g-force may cause tube breakage, with release of sample, droplets, and an aerosol into the centrifuge bowl. Release of these potentially hazardous materials can be avoided by using specially designed sealed containers in which tubes are held during centrifugation. Centrifuge carriers and inserts should be of the size specific to the tubes used. Use of carriers too large or too small for the tube may result in breakage.

Warnings

1. Practice Universal Precautions. Use gloves, gowns, eye protection, other personal protective equipment, and engineering controls to protect from blood splatter, blood leakage, and potential exposure to bloodborne pathogens.
2. Handle all biologic samples and blood collection "sharps" (lancets, needles, luer adaptors and blood collection sets) according to the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury), since they may transmit viral hepatitis, HIV (AIDS), or other infectious diseases. Utilize any built-in used needle protector, if the blood collection device provides one. BD does not recommend resheathing used needles. However, the policies and procedures of your facility may differ and must always be followed.
3. Discard all blood collection "sharps" in biohazard containers approved for their disposal.
4. Transferring a sample collected using a syringe and needle to a tube is not recommended. Additional manipulation of sharps, such as hollow bore needles, increases the potential for needlestick injury.
5. Transferring samples from syringe to an evacuated tube using a non-sharps device should be performed with caution for the following reasons.
 - Depressing the syringe plunger during transfer can create a positive pressure, forcefully displacing the stopper and sample and causing a potential blood exposure.
 - Using a syringe for blood transfer may also cause over or under filling of tubes, resulting in an incorrect blood-to-additive ratio and potentially incorrect analytic results.
 - Evacuated tubes are designed to draw the volume indicated. Filling is complete when vacuum no longer continues to draw, though some tubes may partially fill due to plunger resistance when filled from a syringe. The laboratory should be consulted regarding the use of these samples.
6. If blood is collected through an intravenous (I.V.) line, ensure that line has been cleared of I.V. solution before beginning to fill blood collection tubes. This is critical to avoid erroneous laboratory data from I.V. fluid contamination.
7. Discard blood collection tubes in biohazard containers approved for their disposal.
8. BD Vacutainer® RST Tubes are not recommended for amino acid determinations. Therapeutic drug monitoring (TDM), blood banking and infectious disease performance has not been established, except for anti-CMV IgG and anti-CMV IgM.

ANALYTIC EQUIVALENCY

Evaluations of BD Vacutainer® RST Tubes have been performed for selected analytes on certain instrument platforms. See Table 1. Stability for 24 hours at room temperature has been demonstrated for the listed analytes except LDL (12 hours), progesterone (6 hours), and triglycerides (12 hours). The BD Technical Services Department is available to provide information regarding these studies.

Whenever changing any manufacturer's blood collection tube type, size, handling, processing or storage condition for a particular laboratory assay, the laboratory personnel should review the tube manufacturer's data and their own data to establish/verify the reference range for a specific instrument/reagent system. Based on such information, the laboratory can then decide if changes are appropriate.

REFERENCES

Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS). Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition. Document H1-A5. Wayne, PA: CLSI, 2003.

Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS). Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline – Third Edition. Document H18-A4. Wayne, PA: CLSI, 2009.

TECHNICAL SERVICES

U.S. customers please call BD Global Technical Services at 1.800.631.0174.

Customers outside the U.S. please contact your local BD sales consultant.

All customers: www.bd.com/vacutainer/contact/

TABLE 1**BD Vacutainer® Rapid Serum Tube (RST) Blood Collection Tube****Clinical Evaluations**

Analyte	Instruments
Alanine Aminotransferase.....	7, 8, 9
Albumin	7, 8, 9
Alkaline Phosphatase.....	7, 8, 9
Amylase.....	2, 7, 8
Aspartate Aminotransferase	2, 7, 8
Bilirubin, Direct	2, 7, 8
Bilirubin, Total	7, 8, 9
Blood Urea Nitrogen.....	7, 8, 9
C-Reactive Protein	8, 9
Calcium	7, 8, 9
Carbon Dioxide, Total	2, 7, 8
Chloride	2, 7, 8
Cholesterol.....	6, 7, 8
Complement C3.....	2, 6
Complement C4.....	2, 6
Cortisol.....	3, 10
Creatine Kinase-MB fraction	1, 9
Creatine Kinase, Total	8, 9
Creatinine	7, 8, 9
Cytomegalovirus Antibodies, IgG.....	4, 5
Cytomegalovirus Antibodies, IgM.....	4, 5
Estradiol	3, 10
Ferritin	3, 10
Folate	3, 10
Follicle Stimulating Hormone.....	3, 10
Free Thyroxine.....	3, 10
Free Triiodothyronine	3, 10
Gamma-glutamyltransferase	2, 7, 9
Glucose.....	7, 8, 9
High Density Lipoprotein.....	6, 8, 9
Human Chorionic Gonadotropin.....	3, 10
Immunoglobulin G	2, 6
Immunoglobulin M	2, 6
Iron.....	6, 7, 8
Lactate Dehydrogenase	2, 7, 9
Lipase	2, 7, 8
Low Density Lipoprotein.....	6, 7, 8
Luteinizing Hormone	3, 10
Magnesium	7, 8, 9
Myoglobin	1, 9
Phosphorus	7, 8, 9
Potassium.....	2, 7, 8
Progesterone	3, 10
Rheumatoid Factor.....	2, 6
Sodium.....	2, 7, 8
Testosterone	3, 10
Thyroid Stimulating Hormone.....	1, 3, 10
Total Protein	7, 8, 9
Total Thyroxine.....	1, 3, 10
Total Triiodothyronine	1, 3, 10
Transferrin	2, 7, 9
Triglycerides.....	6, 7, 8
Troponin I	1, 9
Uric Acid.....	2, 7, 9
Vitamin B12	3, 10

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RÉFÉRENCES

Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS). Tubes and additives for venous blood specimen collection; approved standard – fifth edition. Document H1–A5. Wayne, PA: CLSI, 2003.

Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS). Procedures for the handling and processing of blood specimens; approved guideline – third edition. Document H18-A3. Wayne, PA: CLSI, 2004.

SERVICES TECHNIQUES

Veillez contacter le conseiller commercial BD local.
www.bd.com/vacutainer/contact/

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REFERENTIES

Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS). Tubes and additives for venous blood specimen collection; approved standard – fifth edition. Document H1-A5. Wayne, PA: CLSI, 2003.
 Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS). Procedures for the handling and processing of blood specimens; approved guideline – third edition. Document H18-A3. Wayne, PA: CLSI, 2004.

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REFERANSER

Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS). Tubes and additives for venous blood specimen collection; approved standard – fifth edition. Document H1–A5. Wayne, PA: CLSI, 2003.
 Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS). Procedures for the handling and processing of blood specimens; approved guideline – third edition. Document H18–A3. Wayne, PA: CLSI, 2004.

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REF -- Batch Code

LOT -- Catalog Number

STERILE R -- Method of Sterilization Using Irradiation

-- Attention, See Instruction For Use

	EN -- Recyclable
BG --	Може да се подлага на рециклиране
CS --	Recyklovatelné
DA --	Genanvendelig
DE --	Recyclen
EL --	Ανακυκλώσιμο
ES --	Reciclabre
ET --	Ümbertöödeldav.
FI --	Kierrätettävä
FR --	Recyclable
HR --	Može se reciklirati
HU --	Visszaforgatható
IS --	Endurvinnanlegt
IT --	Riciclabile
KK --	Қайтадан өнделіне алынады
KO --	재활용
LT --	Perdirbamas
LV --	Derīgs atbērīzējai pārstrādei
NL --	Recycleerbaar
NO --	Gjenvinnbar
PL --	Podlega utylizacji
PT --	Reciclável
RO --	Reciclabil
RU --	Пригодно для переработки
SK --	Recyklovateľné
SL --	Primerno za reciklažo
SV --	Återanvändbar
TR --	Geri Dönüşümlü
UK --	Підлягає вторинній переробці
ZH --	可回收再利用

EC REP	EN -- Authorized Representative
BG --	Оторизиран представител
CS --	Oprávněný zástupce
DA --	Autoriseret repræsentant
DE --	Autorisierter Vertreter
EL --	Εξουσιοδοτημένος αντιπρόσωπος
ES --	Representante autorizado
ET --	Volitatud esindaja
FI --	Valtuutettu edustaja
FR --	Représentant autorisé
HR --	Ovlašteni predstavnik
HU --	Meghatalmazott képviselő
IS --	Viðurkenndur fulltrúi
IT --	Rappresentante autorizzato
KK --	Ресми өкілі
KO --	인증 대리인
LT --	Įgaliojatis atstovas
LV --	Pilnvarotais pārstāvis
NL --	Erkende vertegenwoordiger
NO --	Autorisert representant
PL --	Autoryzowany przedstawiciel
PT --	Representante autorizado
RO --	Reprezentant autorizat
RU --	Уполномоченный представитель
SK --	Autorizovaný zástupca
SL --	Pooblaščen zastopnik
SV --	Auktoriserad representant
TR --	Yetkili Temsilci
UK --	Вповноважений представник
ZH --	授权代理商

Becton, Dickinson and Company, Franklin Lakes, NJ 07417 USA
 Becton, Dickinson and Company, Bellver Industrial Estate, Plymouth, PL6 7BP UK.
 For reference material, www.bd.com/vacutainer or call 1.800.631.0174 (USA only)
 BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company.
 Also BD Vacutainer® & BD Hemogard™ ©2010 BD. 4015002
 Mfg for BD, Japan 07/2010