

Influenza Quiz

Circle the correct answer

This quiz is an educational tool intended to assist facilities in evaluating their operators' understanding of the QuickVue Influenza Test procedure. This quiz is not intended to be used as sole evidence of operator training or competency. Facilities are responsible for ensuring the quality of the testing performed by their operators. When testing controls or patient specimens, follow the current Package Insert instructions and/or Procedure Card included in the kit.

1. **How often does the manufacturer recommend that external controls be run on the kit?**
 - a) Once every 20 tests
 - b) Once per each new shipment of kits
 - c) Every first sample tested each day
 - d) Once per each untrained operator, once for each new shipment of kits – provided that each different lot in a shipment is tested
2. **At what temperature do the kits need to be stored?**
 - a) 2°C to 8°C (refrigerated)
 - b) 15°C to 30°C (room temperature)
 - c) 2°C to 30°C (refrigerated or room temperature)
 - d) -20°C (frozen)
3. **For proper test performance, which swab should be used when collecting a nasal swab specimen?**
 - a) The swabs provided in the kit
 - b) Any type of swabs
 - c) Any Dacron® swabs
 - d) Cotton-tipped swabs
4. **What is the read time for this kit?**
 - a) 10 minutes
 - b) 15 minutes
 - c) 5 minutes
 - d) 20 minutes
5. **Can the test be read after the read time?**
 - a) No
 - b) Yes, up to 20 minutes after the read time
 - c) Yes, up to 5 minutes after the read time
 - d) Yes, any time after the read time
6. **How long should the swab incubate before adding the dipstick?**
 - a) 0 minutes
 - b) 1 minute
 - c) 10 minutes
 - d) 20 minutes
7. **How long can the swab be stored at room temperature or refrigerated prior to testing?**
 - a) 8 hours
 - b) 12 hours
 - c) 24 hours
 - d) Only a blue control line should be seen
8. **What different sample types are appropriate for testing with this kit?**
 - a) Nasal swab only
 - b) Nasal swab, nasopharyngeal swab, nasal aspirate and nasal wash
 - c) Nasal swab and nasal wash
 - d) Nasal swab, nasal aspirate and nasal wash
9. **How long after adding Extraction Reagent solution into the Extraction Tube is the mixture still good?**
 - a) 8 hours
 - b) 1 minute
 - c) The mixture should be used immediately
 - d) 24 hours
10. **Can transport media be used with this kit?**
 - a) Yes, for all specimen types
 - b) No
 - c) Yes, but only for swab specimens
 - d) Yes, but only for nasal wash specimens

Influenza Quiz

Circle the correct answer

This quiz is an educational tool intended to assist facilities in evaluating their operators' understanding of the QuickVue Influenza Test procedure. This quiz is not intended to be used as sole evidence of operator training or competency. Facilities are responsible for ensuring the quality of the testing performed by their operators. When testing controls or patient specimens, follow the current Package Insert instructions and/or Procedure Card included in the kit.

1. How often does the manufacturer recommend that external controls be run on the kit?
 - a) Once every 20 tests
 - b) Once per each new shipment of kits
 - c) Every first sample tested each day
 - ☒ d) Once per each untrained operator, once for each new shipment of kits – provided that each different lot in a shipment is tested
2. At what temperature do the kits need to be stored?
 - a) 2°C to 8°C (refrigerated)
 - ☒ b) 15°C to 30°C (room temperature)
 - c) 2°C to 30°C (refrigerated or room temperature)
 - d) -20°C (frozen)
3. For proper test performance, which swab should be used when collecting a nasal swab specimen?
 - ☒ a) The swabs provided in the kit
 - b) Any type of swabs
 - c) Any Dacron® swabs
 - d) Cotton-tipped swabs
4. What is the read time for this kit?
 - ☒ a) 10 minutes
 - b) 15 minutes
 - c) 5 minutes
 - d) 20 minutes
5. Can the test be read after the read time?
 - ☒ a) No
 - b) Yes, up to 20 minutes after the read time
 - c) Yes, up to 5 minutes after the read time
 - d) Yes, any time after the read time
6. How long should the swab incubate before adding the dipstick?
 - ☒ a) 0 minutes
 - b) 1 minute
 - c) 10 minutes
 - d) 20 minutes
7. How long can the swab be stored at room temperature or refrigerated prior to testing?
 - ☒ a) 8 hours
 - b) 12 hours
 - c) 24 hours
 - d) Only a blue control line should be seen
8. What different sample types are appropriate for testing with this kit?
 - a) Nasal swab only
 - b) Nasal swab, nasopharyngeal swab, nasal aspirate and nasal wash
 - c) Nasal swab and nasal wash
 - ☒ d) Nasal swab, nasal aspirate and nasal wash
9. How long after adding Extraction Reagent solution into the Extraction Tube is the mixture still good?
 - a) 8 hours
 - b) 1 minute
 - ☒ c) The mixture should be used immediately
 - d) 24 hours
10. Can transport media be used with this kit?
 - ☒ a) Yes, for all specimen types
 - b) No
 - c) Yes, but only for swab specimens
 - d) Yes, but only for nasal wash specimens