

ENZYCLEAN® SUPER CONCENTRATE

Protease Enzyme Low Suds Detergent



Enzyclean® Super Concentrate Protease Enzyme Low Suds Detergent is highly concentrated for use in automatic washing equipment. The protease enzyme in Enzyclean® breaks down organic and inorganic soils, residues, films and deposits and ensures thorough and complete cleaning, leaving instruments moving freely and spot-free.

- **Protease enzyme effectively breaks down proteins.**
- **Super Concentrate formula requires lower dilution, reducing cost.**
- **Compact container makes handling easy and takes up less space.**

Enzyclean® Super Concentrate



**Protease Enzyme
Low Suds Detergent**

Order No.	Size
EBL2.5C	2.5 Gallons (9.5 Liters)

Physical Properties	
Form:	Liquid
Color:	Green
Odor:	Fresh
pH Concentration:	8.7

Uses	Temperature
<input type="checkbox"/> Presoak <input checked="" type="checkbox"/> Automatic Washer	Usage: 60°-140° F (16°-60° C)
<input type="checkbox"/> Manual <input checked="" type="checkbox"/> Automatic Endoscope Reprocessor	Storage: Cool, dry area
<input checked="" type="checkbox"/> Ultrasonic <input type="checkbox"/> Cart Washer	

Product Features	Consumer Benefits
Enzyme:	Protease enzyme effectively breaks down proteins.
Rust Inhibitors:	Reduces likelihood of rust, tarnish or corrosion.
Rinsing Agents:	Rinses spot & residue-free.
Super Concentrated:	Product requires lower dilution, reducing cost.
Compact Container:	Smaller 2.5 gallon container makes handling easy and takes up less space.
Neutral pH:	Product safe and compatible with all metals.
Low Sudsing:	Compatible with most equipment & pump systems.
Bacteriostatic:	If stored properly, not favorable to bacterial growth.

Procedures	Directions for Use
Directions:	Wear appropriate personal protective equipment based on safe work practices.
Ultrasonic Cleaning:	Add 1/20 fl. oz. (1.5 mL) to 1/4 fl. oz. (7.4 mL) per gallon (3.8 L) of water and follow equipment manufacturer's instructions. Recommended temperature range: 90°-140° F (32°-60° C). Lid should be closed during use. Discard solution after each use or when visibly soiled.
Automatic Washing Equipment:	Add 1/20 fl. oz. (1.5 mL) to 1/4 fl. oz. (7.4 mL) per gallon (3.8 L) of water and follow equipment manufacturer's instructions. Recommended temperature range: 90°-140° F (32°-60° C). Validated cleaning is recommended.