

MICROKLEEN™ LS SUPER CONCENTRATE

Low Suds Neutral Detergent



MicroKleen™ LS Super Concentrate Low Suds Neutral Detergent is specially formulated for manual cleaning and for use in ultrasonic washers and automatic washing equipment. It is safe for cleaning glass, plastic, rubber and metal surgical instruments and equipment.

- **Compatible with all equipment & pump systems.**
- **Super concentrate formula requires lower dilution, reducing cost.**
- **Compact container makes handling easy and takes up less space.**

MicroKleen™ LS Super Concentrate

Low Suds Neutral Detergent



Order No.	Size	Physical Properties
DM2.5C	2.5 Gallons (9.5 Liters)	Form: Liquid Color: Blue Odor: Fresh pH Concentration: 8.8

Uses	Temperature
<input type="checkbox"/> Presoak <input checked="" type="checkbox"/> Manual <input checked="" type="checkbox"/> Ultrasonic	Usage: 68°-110° F (20°-43° C) Storage: Cool, dry area
<input checked="" type="checkbox"/> Automatic Washer <input type="checkbox"/> Automatic Endoscope Reprocessor <input checked="" type="checkbox"/> Cart Washer	

Product Features	Consumer Benefits
Rust Inhibitors:	Reduces likelihood of rust, tarnish or corrosion.
Rinsing Agents:	Rinses spot & residue-free.
Neutral pH:	Product safe and compatible with all metals.
Super Concentrated:	Product requires lower dilution, reducing cost.
Compact Container:	Small 2.5 gallon (9.5 L) container makes handling easy and takes up less space.
Low Sudsing:	Compatible with all equipment & pump systems.
Bacteriostatic:	If stored properly, not favorable to bacterial growth.
Environmentally-Friendly:	Biodegradable, phosphate-free, VOC free formula.

Procedures	Directions for Use
Directions:	Wear appropriate personal protective equipment based on safe work practices.
Manual/Ultrasonic Cleaning:	Add 1/20 fl. oz. (1.5 mL) to 1/4 fl. oz. (7.4 mL) per gallon (3.8 L) of warm water. Discard solution after each use.
Automatic Washing Equipment:	Add 1/20 fl. oz. (1.5 mL) to 1/4 fl. oz. (7.4 mL) per gallon (3.8 L) of warm water. Dilutions vary depending on temperature, water quality and equipment. Validated cleaning is recommended.