
# Spinalscopics® Spinal Fluid Cell Count Control / Level 1 & 2

## ENGLISH

### Intended Use

Spinalscopics Spinal Fluid Cell Count Control is intended for monitoring total cell counts performed manually using a hemocytometer to validate quantitation of red and white blood cells in patient CSF samples.

### Summary and Explanation

Microscopic evaluation of cerebrospinal fluid is an important part of the routine analysis of CSF specimens. Along with physical, chemical, immunological, and microbiological examination, total cell count and differential cell count can provide valuable information regarding infection, malignancy and hemorrhage of the central nervous system. Total cell counts are performed in a counting chamber with undiluted CSF.

### Product Description

Spinalscopics Spinal Fluid Cell Count Controls are supplied in two levels, 3 x 3 mL bottles per level per box. The controls are ready-to-use liquid, requiring no reconstitution or dilution. They are prepared in a human protein matrix fortified to target levels with purified chemicals and stabilized human red and white blood cells. Preservatives have been added to inhibit microbial growth.

### Warnings and Precautions

**POTENTIAL BIOHAZARDOUS MATERIAL.** All blood donor units comprising the human source material used in the manufacture of this product have been tested and found non-reactive for Hepatitis B surface antigen, Hepatitis C, and HIV antibodies when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. These samples should be handled according to the Centers for Disease Control's Biosafety Level 2 recommendations.

### Storage and Stability

The Spinalscopics Controls should be stored tightly capped and refrigerated (2°C–8°C) when not in use. Do not freeze. When stored unopened at 2°C–8°C, the controls are stable until the expiration date stated on the label. Once opened, the controls are stable for 6 months when stored at 2°C–8°C between uses. When stored at room temperature (18°C–25°C) the controls are stable for one month even after opening. Discard the controls if there is any evidence of microbial contamination. The level 2 control may appear slightly turbid after mixing.

### Procedure

#### Hemocytometer Cell Count Procedure

1. Remove the control from the refrigerator and replace the cap on the control bottle with the bulb dropper cap included in the control box.
2. Allow the control to come to room temperature (**18°C–25°C**) for at least 15 minutes.
3. Mix the control thoroughly by inverting the bottle **at least 10 times**. Squeeze the bulb in the cap to aspirate and expel the control through the glass dropper attached to the cap (**at least 10 times**) **immediately** prior to use to ensure homogeneity of the content. Thorough mixing with each use is important in order to obtain reproducible results. Avoid foaming.

**Note:** To prevent cells from settling and to ensure homogeneity of sample, the hemocytometer chambers must be charged immediately after thorough mixing. Do not allow cells to settle in bottle prior to charging hemocytometer chambers.

4. Using the glass dropper provided, **charge both sides of the hemocytometer chamber by placing the glass dropper tip alongside cover-slip for each chamber. To prevent overflow chambers, gently squeeze the bulb to slowly charge the chamber until full.**
5. Immediately recap the controls after each use. The controls can be stored for one month at room temperature when not in use. Treat the controls as you would a patient sample.
6. Allow the cells to settle by placing the charged hemocytometer in a humidified chamber for 10 minutes before counting.
7. Count the cells in all nine (9) large squares on both sides of the hemocytometer. **If there are cells on the line, count only those on the top and right side of each square.**
8. Average the number of cells counted on both sides of chamber. Calculate the total number of erythrocytes (RBC) and lymphocytes (WBC) per  $\mu\text{L}$  as follows:

Average cell count (cells/ $\mu\text{L}$ )  $\div 0.9 =$  Total cell count (cells/ $\mu\text{L}$ )

**Note:** Clinical data demonstrates that the majority of leukocytes (WBC) found in Cerebrospinal Spinal Fluid (CSF) are lymphocytes.<sup>1</sup> To more closely match clinical sample, the leukocytes (WBC) population in this product has been lymphocytes enriched. As lymphocytes are similar in size to erythrocytes (RBC), it is important to ensure that laboratory personnel can differentiate lymphocytes (WBC) from erythrocytes (RBC).

**Any deviation from this procedure stated above may introduce variables that can result in values outside specified product insert range. Each laboratory should establish its own precision parameters.**

### Expected Ranges

The expected ranges for red and white blood cells were determined by intralaboratory assays of multiple bottles of the indicated lots of Spinalscopics Control by hemocytometer cell count. The expected ranges are based on the average number of cells per  $\mu\text{L}$  counted in separate assays in the nine large grid squares on both sides of the hemocytometer counting chamber. Use of other systems or protocols may yield differing results. Each laboratory should establish its own precision parameters. A quality control program, Quantrol,® is available from Quantimetrix. Please inquire.

### Limitations

Any future changes made by the manufacturer of a test system may give different values from the indicated range. Detailed information on the limitations of each test system is included in the limitations section of the system manufacturer's package insert. Technical updates can be found on our website. The Quality Control Log will no longer be shipped with the product. You can obtain it from the Quantimetrix website at [www.quantimetrix.com](http://www.quantimetrix.com) or contact Tech Support at (310) 536.0006, option 3.

### References

<sup>1</sup>Burtis C, Ashwood E, ed. Tietz Textbook of Clinical Chemistry, 3rd Ed. W.B. Saunders Company, Philadelphia.

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Quantimetrix Corporation
MDSS GmbH