




## Gelling Fiber Dressing

- High retention capacity to prevent leakage and maceration
- High tensile strength to enable dressing removal in one piece
- Absorbs and retains exudate, blood and bacteria
- Highly absorbent, also under compression therapy
- Soft and conformable which makes it easy to apply

### High retention capacity

- Prevents leakage
- Prevents maceration
- Absorbs and retains exudate, bacteria and blood




### Hydrolock Technology

- Highly absorbent poly-vinyl alcohol fibres
- Excellent tensile strength
- Dressing removal in one piece
- Saves nursing time and increases patient comfort

## How to use Exufiber



1. Cleanse the wound with saline solution or water. Dry the surrounding skin thoroughly.



- 2a. Apply a dry Exufiber dressing to the wound. The dressing should be able to cover the entire wound. It should overlap the dry surrounding skin by at least 1-2 cm for the smaller sizes and 5 cm for larger sizes. The dressing will shrink as it absorbs wound fluid and starts gelling.

- 2b. Loosely pack ribbon or sheet into the wound to allow room for swelling of the dressing.



3. Fixate with an appropriate secondary dressing. Compression therapy may be used in conjunction with Exufiber.

## How Exufiber works

Exufiber is a sterile nonwoven dressing made from highly absorbent polyvinyl alcohol fibers. In contact with wound exudate, Exufiber transforms into a gel that facilitates moist wound healing and ease of removal during dressing change. Exufiber absorbs and retains wound exudate. Exufiber is available both as sheet and ribbon dressings.

## Benefits of Exufiber

- High tensile strength to enable dressing removal in one piece
- High retention capacity to prevent leakage and maceration
- Absorbs and retains exudate, blood and bacteria
- Highly absorbent, also under compression therapy
- Soft and conformable which makes it easy to apply
- Can be left in place for up to 7 days and up to 14 days on donor sites

## Precautions

All wounds should be inspected frequently. In case of signs of clinical infection, consult a health care professional for adequate infection treatment.

Exufiber is not intended for dry wounds or full thickness burns.

If the dressing dries out and is difficult to remove, it should be moistened according to local policies (e.g. with sterile saline or sterile water) and allowed to soak until it lifts easily. It may take several minutes for Exufiber to transform into a gel. Remove the dressing by gently cleansing/flushing

Find out more at [www.molnlycke.co.uk](http://www.molnlycke.co.uk)

Molnlycke Health Care, Arenson Centre, Arenson Way, Dunstable, Bedfordshire, LU5 5UL. Telephone: 0800 917 4918. Email: [info.uk@molnlycke.com](mailto:info.uk@molnlycke.com)

## Intended Use

Exufiber wound dressing is intended to be used on a wide range of exuding wounds:

- Leg and foot ulcers
- Pressure ulcers
- Partial thickness burns
- Surgical wounds
- Donor sites
- Malignant wounds
- Dermal lesions and other external wounds inflicted by trauma

Configurations available both for flat and cavity wounds



Exufiber assortment (Sterile packed)

Art. no	Size cm	Pieces Per Inner	PIP Code	NHS Code
603300	5x5	10	394-4337	TBC
603301	10x10	10	394-4345	TBC
603302	15x15	10	394-4352	TBC
603305	4,5x10	10	394-4402	TBC
603306	4,5x20	10	394-4386	TBC
603307	4,5x30	10	394-4378	TBC
603304	2x50	5	394-4394	TBC
603303	20x30	5	394-4360	TBC