




## LIGHTWEIGHT, TUBULAR DRESSING RETENTION

New Tubifast 2-Way Stretch is the world's first and only lightweight tubular bandage with both radial and longitudinal stretch.<sup>1</sup> This unique mode of action allows patients complete freedom of movement, with added comfort. Tubifast is ideal for dressing retention or as a skin covering for any part of the body. It can also be used for patch wrapping and as an under-cast stockinette. Tubifast 2-Way Stretch holds dressings securely in place without constriction or compression. It is quick and easy to use, with less risk of fraying<sup>2</sup> than traditional tubular bandages.



## How Tubifast 2-Way Stretch works

New Tubifast 2-Way Stretch is manufactured from viscose with very fine elastane threads knitted into the fabric radially and longitudinally, to provide light elasticity.

Tubifast 2-Way Stretch holds dressings securely, without constriction or compression. No pins or tapes are necessary. Its two-way stretch allows patients complete freedom of movement and added comfort.


Tubifast 2-Way Stretch is available in a range of sizes designed to use on any part of the body. Each size is identified by a coloured line knitted into the bandage. It's quick and easy to use: cut to size and stretch over the dressing for an even, non-constrictive fit.

## Areas of use

Tubifast 2-Way Stretch can be used as a dressing retention and skin covering for any part of the body. It can also be used for patch wrapping and as an under-cast stockinette.

Because of its tubular construction, Tubifast 2-Way Stretch is particularly suitable for holding dressings in place on difficult areas such as stumps of amputee patients.


## How to use Tubifast 2-Way Stretch



1. Choose the correct size using the quick reference colour coding, and cut to length.



2. Stretch Tubifast 2-Way Stretch and ease gently into place over the affected area.



3. Position Tubifast 2-Way Stretch to cover the dressing. Illustrated, step-by-step instructions are available for more advanced applications.

## Benefits of Tubifast 2-Way Stretch

- Tubifast 2-Way Stretch holds dressings securely, without constriction or compression.
- Its radial and longitudinal stretch allows patients complete freedom of movement.
- A variety of lengths are available, ensuring cost-effectiveness by minimizing waste.
- Quick and easy to use: simply cut to size and ease gently into place over the dressing for an even, non-constrictive fit.
- Available in a range of quick reference, colour-coded sizes to fit everything from small limbs to adult trunks.
- It has less fray so there is less risk of fibres in the wound<sup>2</sup>.

## Fibre content

Viscose, Elastane and Polyamide

## Washing Instructions



### Tubifast 2-Way Stretch Assortment

Art. No.	Product	Width x Length
<b>RED LINE</b>		
2434	Small limbs	3.5cm x 10m
<b>GREEN LINE</b>		
2436	Small and medium limbs	5cm x 10m
<b>BLUE LINE</b>		
2438	Large limbs	7.5cm x 10m
<b>YELLOW LINE</b>		
2440	Extra-large limbs, heads, children's trunks	10.75cm x 10m
<b>PURPLE LINE</b>		
2444	Adult trunks	20cm x 10m

Tubifast is packaged in single packs

Natural rubber latex is not a constituent of Tubifast 2-Way Stretch or its packaging.

<sup>1</sup> Published data shows that Tubifast 2-Way Stretch exceeds the minimum requirement of longitudinal stretch for the actions of sitting, bending and flexing. Physical testing of textiles, B. P. Saville, Woodhead Publishing, 1999.

<sup>2</sup> Data on file.

Mölnlycke Health Care  
1-800-494-5134

The Mölnlycke Health Care name and logo and Tubifast 2-Way Stretch® are registered trademarks of Mölnlycke Health Care AB.



CA171500807EN