

■ Collection and dispensing from a single plastic tube!

■ Fill easily by capillary action to preset volume. No aspiration, air bubbles, or wasted samples.

■ Dispense easily by pinching the unfilled portion of the tube. No bulbs, dispensers or other mechanical devices needed. Ideal for fingerstick sampling and dispensing.



Air vent regulates air flow and volume

MICRO SAFE[®] TUBE

One-piece collection and dispensing

medicore

Ordering information: Medicore[™]

Reorder Number: MS 1015

www.MediCore.com

800-327-8894

Nashville, TN

Printed in U.S.A.


All Rights Reserved.

© 2002 Safe-Tec Clinical Products, LLC

U.S. Patent No. 5,460,782


MICROSAFE® Tube

An air vent is strategically positioned on the sidewall of the tube to provide automatic air venting and sample volume control.




CAUTION!

Filling is automatic: Never squeeze the tube while sampling.



Step 1: Hold the tube horizontally, and touch the tip of the MICROSAFE Tube to the blood sample.

Capillary action will automatically draw the sample to the fill line and stop.



Step 2: To expel the sample, align the tip of the tube with the sample target and squeeze the bulb.

If a sample won't expel, you probably didn't allow the tube to fill all the way. Touch the tip to the blood sample again and allow it to fill completely. Then align the tip with the sample target and squeeze the bulb.