


SORBACT® FOAM DRESSING

Infection preventive foam dressing



REDUCES BIOBURDEN IN MODERATELY EXUDING WOUNDS

Sorbact® Foam Dressing is a sterile, bacteria and fungi binding foam dressing. It consists of a green Sorbact® wound contact layer combined with a polyurethane foam with a semi-permeable backing. Sorbact® Foam Dressing absorbs and retains exudate, thereby reducing the risk of maceration and enabling a moist wound environment. The semi-permeable film backing allows excess fluid to evaporate.

Product	Ref. no.	Size	Pcs/pack
Sorbact® Foam Dressing	98310	10x10 cm	10
Sorbact® Foam Dressing	98315	15x15 cm	10
Sorbact® Foam Dressing	98320	10x20 cm	10




INTENDED USE

Sorbact® Foam Dressing is intended for use in management of clean, colonized, contaminated or infected wounds with low to moderate exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

INSTRUCTIONS FOR USE

1. Prepare the wound according to local clinical practice.
2. Select an appropriate dressing size for the wound. The wound pad should overlap the wound margins by at least 2 cm.
3. Remove the dressing from the pouch using an aseptic technique.
4. If cutting the dressing, use an aseptic technique. Discard any open and unused dressing.
5. Apply the dressing. Ensure that the green side of the dressing comes into direct contact with the complete wound surface.
6. Apply an appropriate dressing fixation.
7. The dressing change frequency depends on exudate levels and overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.



SORBACT® – INFECTION PREVENTIVE DRESSINGS

Sorbact® dressings employ a natural mechanism to bind microorganisms to their surface. The microorganisms are removed when the dressing is changed and the bioburden is thereby decreased.