US

CA

AU

CIDEZYME® XTRA

Multi Enzymatic Detergent for Automated Systems and Manual Cleaning

INSTRUCTIONS FOR USE

CIDEZYME® XTRA Multi Enzymatic Detergent for Automated Systems and Manual Cleaning when used as directed is a low-foaming, mild pH cleaning solution which combines chemical and multi-enzymatic (Cellulase, Lipase, Protease and Amylase) actions to rapidly remove medical soils (e.g. blood). The product also contains corrosion inhibitors to ensure compatibility with medical devices. CIDEZYME XTRA Multi Enzymatic Detergent is uniquely formulated for use in automatic cleaning systems (e.g. endoscope cleaner and reproprocessors, automatic endoscope reproprocessors) as well as for manual cleaning.

**RECOMMENDED CONCENTRATIONS**  
For use in automatic cleaning systems, minimum 4 ml per 1 liter (1/2 oz per gallon) of water or per automatic cleaning system manufacturer’s instructions.  
For manual cleaning, minimum 8 ml per 1 liter (1 oz per gallon) of water.

**CAUTION**  
As with all chemicals, good industrial hygiene practice should be followed. For further information refer to the Material Safety Data Sheet.

**PLEASE NOTE**  
• Confirm CIDEZYME XTRA Multi Enzymatic Detergent compatibility with the instruments and equipment prior to use.  
• Do not add bleach or other chemicals.  
• Soak and clean instrument immediately after use, avoid allowing the medical soil to become dried.  
• Minimum in-use temperature for detergent is 15°C.  
• Unless the automatic cleaning system is specifically designed for elimination of manual cleaning, it will be necessary to pre-clean instruments with brushes or soft cloths as recommended by the instrument manufacturer.  
• CIDEZYME XTRA Multi Enzymatic Detergent is not intended as a solution for holding or storing equipment.

REF 22591

CIDEZYME® XTRA

Multi Enzymatic Detergent for Automated Systems and Manual Cleaning

ADVANCED STERILIZATION PRODUCTS

Division of Ethicon, Inc.

a  company

33 Technology Drive, Irvine, CA 92618

US

Irvine, CA 92618

BR

Rodovia Presidente Dutra, Km 154  
S.J. Campos-S.P 12240-908

CA

Johnson & Johnson Medical Products  
Markham, ON, L3R 0T5

RU

Johnson & Johnson Medical ASP LLC  
121614 Moscow, Russia

AU

Johnson & Johnson Medical Pty Ltd  
1-5 Khartout Road  
North Ryde, NSW 2113

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**DIRECTIONS FOR USE**  
**Automatic Cleaning & Disinfecting System (e.g. automatic endoscope reproprocessors or endoscope cleaner and reproprocessors)**  
1) Connect bottle of CIDEZYME XTRA Multi Enzymatic Detergent to the machine or dispense into detergent reservoir as per system manufacturer’s instructions.  
2) Select a rinse cycle on an automatic system that has been validated for use with this product.  
3) Ensure that the automatic rinse cycle selected will thoroughly rinse the medical device including all lumens with large volume of water equivalent to the reusable device manufacturer’s recommendations.

**Manual Cleaning**  
1) Attach the pump to the bottle and secure tightly.  
2) Add minimum 8 ml of CIDEZYME XTRA Multi Enzymatic Detergent liquid concentrate per 1 liter of water (1 oz per gallon).  
3) Soak instruments for a minimum of 3 minutes. Extend soak time or use warm water for equipment with dried-on organic matter. Flush solution through all channels. Soaking for long periods of time may cause discoloration of anodized aluminum.  
4) Thoroughly rinse the equipment, aspirating water through all channels to remove traces of CIDEZYME XTRA Multi Enzymatic Detergent.  
5) Roughly dry equipment, including channels. Equipment is now ready for reprocessing.  
6) Discard used solution after each use.

**STORAGE**  
CIDEZYME XTRA Multi Enzymatic Detergent should be stored in its original container at controlled temperature 15°- 30°C (59°- 86°F). Replace cap when not in use.

**DISPOSAL**  
In-use concentration of CIDEZYME XTRA Multi Enzymatic Detergent can be disposed directly to drain after use or per hospital policy. After rinsing with water, the empty containers can be disposed of as non-hazardous waste or per hospital policy. Use by expiration date on the immediate container.

**RISK/SAFETY INFORMATION**  
**S24/25** - Avoid contact with skin and eyes.  
Wear appropriate personal protective equipment.  
Skin: Prolonged or repeated contact may cause drying of skin with risk of dermatitis (red, cracked skin). Symptoms are generally alleviated when exposure ends.  
Eyes: May cause light conjunctivitis (redness), stinging or excess tearing. Do not ingest: may cause harm.

Contains Cellulase, Lipase, Protease and Amylase. May produce an allergic reaction.

HOW SUPPLIED	
<b>Reorder</b>	<b>Description</b>
22591	2 X 3.8L (1 gallon)

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