


**OSTEOMED**

A TRADITION OF QUALITY AND INNOVATION

# NEUROPACK<sup>TM</sup>

STERILE NEURO FIXATION SYSTEM
























- STERILE TITANIUM PLATES & SCREW KITS
- THINNEST PROFILE AVAILABLE 0.25mm
- COMPLIMENTS EXISTING FAST-FLAP<sup>TM</sup> SYSTEM
- MEETS TRACEABILITY NEEDS
- ALL KITS INCLUDE THE PATENTED AUTO-DRIVE<sup>®</sup> SELF DRILLING SCREWS WHICH ELIMINATE THE NEED FOR PRE-DRILLING

PRODUCTS THAT IMPROVE PEOPLE'S LIVES

# NEUROPACK™

## STERILE NEURO FIXATION SYSTEM

part number	content	qty	description
218-0070-SP		2	Medium 2 Hole Straight Plate
		1	Long 2 Hole Straight Plate
		6	1.6mm x 4mm Auto-Drive® Screw
218-0072-SP		2	Medium 2 Hole Straight Plate
		1	Long 2 Hole Straight Plate
		6	1.6mm x 4mm Auto-Drive® Screw
		1	Manual Screwdriver Assembly
218-0074-SP		3	Medium 2 Hole Straight Plate
		6	1.6mm x 4mm Auto-Drive® Screw
218-0076-SP		2	Medium 2 Hole Straight Plate
		4	1.6mm x 4mm Auto-Drive® Screw
218-0078-SP		1	Large Bur Hole Plate
		6	1.6mm x 4mm Auto-Drive® Screw
218-0174-SP		3	1.6mm, Low Profile Medium 2 Hole Straight Plate
		6	1.6mm x 4mm Auto-Drive® Screw
218-0176-SP		2	1.6mm, Low Profile Medium 2 Hole Straight Plate
		4	1.6mm x 4mm Auto-Drive® Screw
218-0178-SP		1	1.6mm, Low Profile Large Bur Hole Plate
		6	1.6mm x 4mm Auto-Drive® Screw
218-0075-SP		3	Profile 0, 2 Hole Straight Plate
		6	1.6mm x 4mm, Profile 0 Auto-Drive Screw



3885 Arapaho Road  
Addison, Texas 75001  
Customer Service: 800.456.7779  
Fax: 800.390.2620

www.osteomed.com

P/N 030-1582 Rev. B



A COLSON ASSOCIATE

Refer to **Product Insert** for specific indications and **Instructions For Use**.