



# PROFILE 0

LOW PROFILE NEURO PLATING SYSTEM

RETHINKING  
POSSIBILITIES  
RESHAPING  
LIVES

At OsteoMed, we work alongside surgeons who devote their lives to improving those of others. Our goal: rethinking the possibilities that have the power to reshape lives.

With a spirit of innovation and a relentless focus on quality, we offer surgeons a wide range of choices for neurosurgical and craniomaxillofacial procedures.









Patented Auto-Drive screw with friction fit technology

# PROFILE

LOW PROFILE NEURO PLATING SYSTEM




*Rethinking Possibilities, Reshaping Lives*

# Lowest Profile

## Neurosurgical Plating System

Profile Zero is an innovative addition to the OsteoMed Neurosurgical plating product offering. At .25mm plate profile, the system virtually eliminates patient palpability. The uniquely engineered plate and screw construct offers strength for cranial fixation with the lowest profile on the market. With a variety of plate designs, Profile Zero provides the neurosurgeon with many options in the OR. Additionally, the updated tray design reduces the system's "footprint" in the surgical theater. Coupled with the plates and screws is a new ergonomic screwdriver handle which is functional and aesthetically appealing.



- 
- .25mm plate-to-screw profile
  - Beveled plate edges for reduced palpability
  - Plates easily contourable to patient anatomy
  - Customizable VisiDisk
  - Reduced OR footprint



# Profile 0 Implants & Meshes



218-0201 Profile Zero 2 Hole Plate



218-0203 Profile Zero 4 Hole Straight Plate Long



218-0205 Profile Zero Y Plate



218-0207 Profile Zero 9mm Square Plate



218-0209 Profile Zero 2 X 4 Hole Quad Plate




218-0210 Profile Zero Burr Hole Plate Small (Ø14.5mm)



218-0211 Profile Zero Large Burr Hole Plate (Ø22mm)



218-0212 Profile Zero X-Large Burr Hole Plate (Ø26mm)



218-0285 Profile Zero 85mm X 50mm X .25mm Mesh



218-0202 Profile Zero 4 Hole Straight Plate



218-0204 Profile Zero 16 Hole Straight Plate



218-0206 Profile Zero Double Y Plate



218-0208 Profile Zero 2 X 3 Hole Quad Plate



218-0213 Profile Zero Shunt Plate, Groove




218-0214 Profile Zero 25mm Gap Plate



218-0215 Profile Zero 35mm Gap Plate



218-0216 Profile Zero 45mm Gap Plate



218-0290 Profile Zero 65mm X 90mm X .25mm Mesh

# Instruments




220-0262 Blue Metal Screwdriver Body



220-0263 Profile Zero 1.6mm Screwdriver Shaft, Manual



220-0260 Small Grasping Forceps



220-0261 Mesh Cutter, Angled



220-0265 Visidisk, Profile Zero



222-0004 1.1mm Pilot Drill, 4mm Stop, J-Latch



222-0006 1.1mm Pilot Drill, 6mm Stop, J-Latch



450-0600 OsteoDriver™ Disposable AutoDriver™

# Screws



218-1604 1.6mm X 4mm Profile Zero AutoDrive® Screw



218-1605 1.6mm X 5mm Profile Zero Autodrive® Screw



218-1904 1.9mm X 4mm Low Profile Emergency Screw



## OSTEOMED

3885 Arapaho Rd.  
Addison, TX 75001  
Customer Service: 800.456.7779  
Customer Service Fax: 800.390.2620

[www.osteomed.com](http://www.osteomed.com)

P/N 030-1731 Rev.C

