



# smartflex

PEDIATRIC CRANIAL EXPANDER



*Rethinking Possibilities, Reshaping Lives*

# smartflex

PEDIATRIC CRANIAL EXPANDER

## MANUAL CRANIAL SPRING



### SmartFlex Cranial Spring

- Pre-formed based on an extensive analysis of infant skulls
- Medical grade stainless steel
- Adjustable to accommodate various skull anatomies
- Beveled foot plate reduces the possibility of dural tear
- Foot plates are designed to fit securely under the skull bone



## SMARTFLEX SPRING CLIP



### SmartFlex Spring Clip

- Selection of 10 newton forces to accommodate patient anatomy
- Low profile, ergonomic design to facilitate a minimally invasive surgical approach.
- Spring Sleeves secures the springs during implantation.
- Push-button attachment secures the spring clip to the spring inserter instrument.

## SMARTFLEX INSERTER




### SmartFlex Inserter

- Push-Button Feature secures the spring clip to the spring inserter instrument
- Ergonomic Handle
- For comfort and ease of use spring loaded trigger





*Rethinking Possibilities, Reshaping Lives*



 A COLSON ASSOCIATE

**OSTEOMED**  
3885 Arapaho Rd.  
Addison, TX 75001  
Customer Service: 800.456.7779  
Customer Service Fax: 800.390.2620

[www.osteomed.com](http://www.osteomed.com)  
030-1818 Rev.A