

OSTEOMED

A TRADITION OF QUALITY AND INNOVATION

EXTREMIFIX™

Cannulated Screw System



Cannulated




Lesser Digital
Arthrodesis



Headless

Superior Bone Purchase
Excellent Compression
Intra-operative Choices
Four Systems in One




PRODUCTS THAT IMPROVE PEOPLE'S LIVES

ExtremiFix Cannulated Screw System is the only system that offers 2.0, 2.4, 3.0 and 4.0mm cannulated lag screws in both headless and headed screw forms to treat a variety of indications in the upper and lower extremity.



Reverse cutting flute to facilitate in screw extraction




Titanium, Self-drilling & Self-tapping screws provide superior compression



One of the lowest profile cannulated screw heads on the market which helps to minimize soft-tissue irritation

ExtremiFix Cannulated Screw System is designed with the surgeon in mind providing multiple intra-operative choices. OsteoMed's ExtremiFix Headless and Cannulated Screws deliver the same quality and innovation you have come to expect.




Bone Clamp w/ K-Wire Guide:

Fracture stabilization and insertion of a k-wire across the fracture line is easily made with this new instrument



Tissue Elevator:

Use to perform soft tissue release and minimizes the amount of soft tissue dissection and trauma to a particular area



Quick Release:

A new instrument to help aid in changing surgically powered drills, drivers, and countersinks in seconds



Drill Guide:

New double sided drill guide assembly for both 2.0/2.4mm and 3.0/4.0mm cannulated drills for precision drilling and the protection of soft tissue

EXTREMIFIX™

Refer to **Product Insert** for specific indications and **Instructions For Use**.




Instruments and Accessories

PART NUMBER DESCRIPTION

316-1400	ExtremiFix™ Tray
316-0035	2.0/2.4/3.0/4.0mm Cannulated Drill Guide
316-0036	4"/6" Cannulated Depth Gauge, Double Sided
316-0038	Cannulated Screwdriver Body, Swivel
316-0040	Cannulated Screwdriver Body, Ratchet
320-0041	Cannulated Screwdriver Body, Static
316-0039	2.0/2.4/3.0/4.0mm Cannulated Screw Remover
316-0042	Screw Forceps
316-0043	11mm Tissue Elevator
316-0044	13mm Tissue Elevator
316-0047	15mm Tissue Elevator
316-0045	Bone Clamp w/K-wire Guide
316-0046	Bone Clamp Curved Lobster Claw

Quick Release Instruments:

320-0417	Quick Release Adapter
316-0013	2.0/2.4mm Cannulated Drill, Quick Release
316-0014	2.0/2.4mm Cannulated Drill, Long, Quick Release
316-0015	3.0/4.0mm Cannulated Drill, Quick Release
316-0211	2.0/2.4mm Cannulated Countersink, Quick Release
316-0201	3.0/4.0mm Cannulated Countersink, Quick Release
316-0132	2.0/2.4mm Cannulated Driver Stem, Quick Release
316-0232	3.0/4.0mm Cannulated Driver Stem, Quick Release
316-0233	3.0/4.0mm Cannulated Driver Stem, Short, Quick Release
316-0304	2.0/2.4mm Headless Screw, Tri-Lobe Driver, Short, Quick Release
316-0305	2.0/2.4mm Headless Screw, Tri-Lobe Driver, Long, Quick Release
316-0314	3.0/4.0mm Headless Screw, Tri-Lobe Driver, Short, Quick Release
316-0315	3.0/4.0mm Headless Screw, Tri-Lobe Driver, Long, Quick Release
316-0320	2.0/2.4mm Proximal Cortex Drill, Quick Release
316-0321	3.0/4.0mm Proximal Cortex Drill, Quick Release



OSTEOMED

3885 Arapaho Road
Addison, Texas 75001
Customer Service: 800.456.7779
Main: 972.677.4600
Fax: 800.390.2620

www.osteomed.com

P/N 030-1611 Rev.A



A COLSON ASSOCIATE