

OSTEOMED

A TRADITION OF QUALITY AND INNOVATION

CANNULATED SCREW SYSTEM™



CANNULATED SELF-DRILLING LAG SCREWS

- PREMIUM CANNULATED SYSTEM FOR THE PODIATRIST, ORTHOPEDIC FOOT AND ANKLE SPECIALIST & ORTHOPEDIC HAND SURGEON
- INNOVATIVE SELF-DRILLING AND SELF-TAPPING TITANIUM¹ LAG SCREWS ENSURE SUPERIOR COMPRESSION WHILE REDUCING OPERATIVE TIME
- LOW PROFILE HEAD PREVENTS SOFT TISSUE IRRITATION
- COMPREHENSIVE SYSTEMS REQUIRE MINIMAL INSTRUMENTATION

PRODUCTS THAT IMPROVE PEOPLE'S LIVES

1. TITANIUM PER ASTM F136

CANNULATED SCREW S

From the leaders in Extremity fixation... comes a Cannulated System for the quality and innovation you have come to expect with the lowest profile head on the market. Synonymous with OsteoMed. These self-drilling and self-tapping lag screws ensure that the instrument design assures ease of use.

316-0106 .035 x 6" Guide Wire
316-0107 .035 x 4" Guide Wire
316-0003 .045 x 4" Guide Wire

316-0004 .045 x 6" Guide Wire
316-0006 .045 x 6" Partially Threaded Guide Wire

316-0119 .035 x 4" Guide Wire - Double Trocar
316-0004 .045 x 6" Guide Wire - Double Trocar

316-0118 .062 x 6" Guide Wire - Double Trocar

316-0116 2.0/2.4mm Cannulated Drill - Extra Long

316-0105 2.0/2.4mm Cannulated Drill
316-0005 3.0/4.0mm Cannulated Drill

316-0102 2.0/2.4mm Cannulated Screwdriver Shaft
316-0002 3.0/4.0mm Cannulated Screwdriver Shaft - Long
316-0007 3.0/4.0mm Cannulated Screwdriver Shaft - Short

316-0109 2.0/2.4mm Screw Remover
320-4040 3.0/4.0mm Screw Remover

316-0000 Cannulated Depth Gauge - 6"
316-0108 Cannulated Depth Gauge - 4"

316-0101 2.0/2.4mm Cannulated Screw Countersink
316-0001 3.0/4.0mm Cannulated Screw Countersink

316-0226 Washer, 4.0mm Cannulated Screws
316-0225 Washer, 3.0mm Cannulated Screws
316-0126 Washer, 2.4mm Cannulated Screws
316-0125 Washer, 2.0mm Cannulated Screws



316-0103
2.0/2.4mm
TAPERLOCK™
Driver body
with ratchet



320-0102
Bone Clamp


INSTRUMENTATION FOR 2.0MM 2.4MM 3.0MM 4.0MM

319-20XX 2.0mm Cannulated Lag Screw Series, 6mm-42mm
319-24XX 2.4mm Cannulated Lag Screw Series, 6mm-50mm
319-30XX 3.0mm Cannulated Lag Screw Series, 10mm-40mm
319-40XX 4.0mm Cannulated Lag Screw Series, 12mm-52mm



the discriminating surgeon. OsteoMed's Cannulated Screws deliver the same
the market and effective, self-drilling, self-tapping capabilities that are synon-
e superior compression while reducing operative time. Surgeon driven instru-

INSTRUMENTATION FOR 6.5mm 7.3mm



319-6XXX-SP 6.5mm Cannulated Screw Series
319-7XXX-SP 7.3mm Cannulated Screw Series
Screw Specifications: (5mm increments)
6.5mm x 40mm-120mm, 20mm thread length
6.5mm x 45mm-120mm, 32mm thread length
7.3mm x 40mm-120mm, 20mm thread length
7.3mm x 45mm-120mm, 32mm thread length

CANNULATED SCREW SYSTEM™

Refer to **Product Insert** for specific indications and **Instructions For Use**.



316-0110
2.0/2.4mm Cannulated
Screw Organizer Block



316-0033 6.5mm / 7.3mm Cannulated Tray

316-0010
3.0/4.0mm Cannulated
Screw Organizer Block



CANNULATED SELF-DRILLING LAG SCREWS

2.0mm • 2.4mm • 3.0mm • 4.0mm

316-0106	.035 x 6" Guide Wire
316-0107	.035 x 4" Guide Wire
316-0003	.045 x 4" Guide Wire
316-0004	.045 x 6" Guide Wire
316-0006	.045 x 6" Partially Threaded Guide Wire
316-0105	2.0/2.4mm Cannulated Drill
316-0005	3.0/4.0mm Cannulated Drill
316-0102	2.0/2.4mm Cannulated Screwdriver Shaft
316-0002	3.0/4.0mm Cannulated Screwdriver Shaft - Long
316-0007	3.0/4.0mm Cannulated Screwdriver Shaft - Short
316-0109	2.0/2.4mm Screw Remover
320-4040	3.0/4.0mm Screw Remover
316-0101	2.0/2.4mm Cannulated Screw Countersink
316-0001	3.0/4.0mm Cannulated Screw Countersink
316-0000	Cannulated Depth Gauge - 6"
316-0108	Cannulated Depth Gauge - 4"
316-0010	3.0/4.0mm Cannulated Screw Organizer Block
320-0102	Bone Clamp
316-0103	2.0/2.4mm TAPERLOCK™ Driver body with ratchet
316-0125	Washer, 2.0mm Cannulated Screws
316-0126	Washer, 2.4mm Cannulated Screws
316-0225	Washer, 3.0mm Cannulated Screws
316-0226	Washer, 4.0mm Cannulated Screws
319-20XX	2.0mm Cannulated Lag Screw Series, 6mm-42mm
319-24XX	2.4mm Cannulated Lag Screw Series, 6mm-50mm
319-30XX	3.0mm Cannulated Lag Screw Series, 10mm-40mm
319-40XX	4.0mm Cannulated Lag Screw Series, 12mm-52mm

CANNULATED SELF-DRILLING LAG SCREWS

6.5mm • 7.3mm

316-0023	Partially Threaded Guide Wire, 2.8mm X 300mm
316-0019	Cannulated Tap for 6.5mm Screws
316-0031	Cannulated Tap for 7.3mm Screws
316-0016	Cannulated Drill, 6.5mm / 7.3mm
316-0017	Depth Gauge, 6.5mm / 7.3mm
316-0018	Adjustable Parallel Drill Guide, 6.5mm / 7.3mm
316-0021	Cannulated Non-slip Screwdriver Body, 6.5mm / 7.3mm
316-0022	Cannulated Screwdriver Shaft, 6.5mm / 7.3mm
316-0024	Screw Extractor
316-0025	Washer for 6.5mm / 7.3mm Screws
316-0032	Ratchet Cannulated Screwdriver Body, 6.5mm / 7.3mm
316-0020	Cannulated Countersink, 6.5mm / 7.3mm
316-0026	Trocar, 2.8mm
316-0027	Guide Pin Sleeve, 6.5mm / 7.3mm
316-0028	Drill Sleeve, 6.5mm / 7.3mm
316-0029	Screw Sleeve, 6.5mm / 7.3mm
316-0030	Washer Sleeve, 6.5mm / 7.3mm
316-0033	6.5mm / 7.3mm Cannulated Tray
319-6XXX-SP	6.5mm Cannulated Screw Series, 40mm-120mm 20mm or 32mm thread length
319-7XXX-SP	7.3mm Cannulated Screw Series, 40mm-120mm 20mm or 32mm thread length



3885 Arapaho Road
Addison, Texas 75001
Customer Service: 800.456.7779
Main: 972.677.4600
Fax: 800.390.2620
www.osteomedcorp.com

P/N 030-1437 Rev.C



A COLSON ASSOCIATE