

enCompass


Lessers Metatarsal Resurfacing Implant



Rethinking Possibilities, Reshaping Lives

established **DESIGN** **expands** to the **lesser** Metatarsals

- Low profile one-piece implant promotes minimal bone resection
- Four finned stem and reverse cutting teeth prevent implant rotation and pullout
- Titanium plasma spray and hydroxyapatite coating encourages osseointegration
- Streamlined technique and calibrated instrumentation allows for precise decompression of the joint



OSTEOMED

3885 Arapaho Rd.
Addison, TX 75001

Customer Service: 800.456.7779

Customer Service Fax: 800.390.2620

www.osteomed.com

P/N 030-1767 Rev.B



A COLSON ASSOCIATE