

OSTEOMED

A TRADITION OF QUALITY AND INNOVATION

MINCRO™

BONE GRAFT REGENERATION FIXATION KIT



SIMPLE FIXATION, IN A COMPACT KIT

- COMPACT AND COMPLETE KIT FOR ONLAY BONE GRAFTING AND GUIDED TISSUE REGENERATION PROCEDURES
- COST-EFFECTIVE SYSTEM THAT CAN BE CUSTOMIZED TO YOUR PRACTICE
- TITANIUM ALLOY SCREWS PROVIDE SUPERIOR BONE PURCHASE AND STRENGTH
- 1.6MM AUTO-DRIVE™ SELF-DRILLING SCREWS ARE AVAILABLE FOR MEMBRANE TACKING AND ONLAY BONE GRAFTING
- 1.2MM AUTO-DRIVE™ SELF-DRILLING SCREWS ARE AVAILABLE FOR MEMBRANE TENTING, MEMBRANE TACKING, AND ONLAY BONE GRAFTING


PRODUCTS THAT IMPROVE PEOPLE'S LIVES

MINCRO™

BONE GRAFT REGENERATION FIXATION KIT

ONLAY BONE GRAFT AND GUIDED TISSUE REGENERATION FIXATION KIT

Refer to [Product Insert](#) for specific indications and [Instructions For Use](#).




1.2mm Membrane Tenting Screws

Large surface area head effectively supports membrane
Maintains volume in the graft site

1.2mm and 1.6mm Membrane Tacking Screws

Securely stabilizes membrane
Low profile head minimizes soft-tissue perforation



1.6mm Bone Graft Fixation Screws

Exceptional bone purchase
Easy removal without risk of stripping or breaking

GUIDED TISSUE REGENERATION

1.2mm Membrane Tenting Screws

205-1206	1.2mm x 6mm Tenting Screw
205-1208	1.2mm x 8mm Tenting Screw
205-1210	1.2mm x 10mm Tenting Screw

1.2mm Standard Membrane Tacking Screw

204-1202	1.2mm x 2mm Screw
204-1203	1.2mm x 3mm Screw
204-1204	1.2mm x 4mm Screw
204-1206	1.2mm x 6mm Screw

1.6mm Standard Membrane Tacking Screw

208-1604	1.6mm x 4mm Screw
208-1605	1.6mm x 5mm Screw
208-1606	1.6mm x 6mm Screw

1.6mm AUTO-DRIVE™ Membrane Tacking Screw

211-1604	1.6mm x 4mm AUTO-DRIVE™ Screw
211-1606	1.6mm x 6mm AUTO-DRIVE™ Screw

Plates

212-0034	Micro Straight Plate
212-0041	Short 2x4 Hole Quad Plate
220-0021	2.0/2.4mm Screwdriver Shaft, Manual

INSTRUMENTATION

220-0060	Mincro Screw Organizer Block
220-0019	TaperLock™ Screwdriver Body
220-0069	Membrane Punch
220-0061	1.6mm Screwdriver Shaft, Manual
220-0067	1.6mm Screwdriver Shaft, Dental Latch
220-0064	1.3mm Pilot Drill, Dental Latch
220-0022	1.2mm Screwdriver Shaft, Manual
220-0068	1.2mm Screwdriver Shaft, Dental Latch
220-0066	1.0mm Pilot Drill, 4mm Stop, Dental Latch
220-0065	1.0mm Pilot Drill, 8mm Stop, Dental Latch
220-0165	1.0mm Pilot Drill, 12mm Stop, Dental Latch

BONE GRAFT FIXATION

1.2mm Standard Bone Graft Fixation Screws

204-1206	1.2mm x 6mm Screw
204-1208	1.2mm x 8mm Screw
204-1210	1.2mm x 10mm Screw
204-1212	1.2mm x 12mm Screw
204-1214	1.2mm x 14mm Screw

1.6mm AUTO-DRIVE™ Bone Graft Fixation Screws

211-1606	1.6mm x 6mm AUTO-DRIVE™ Screw
211-1608	1.6mm x 8mm AUTO-DRIVE™ Screw
*201-1610	1.6mm x 10mm Partially Threaded AUTO-DRIVE™ Screw
*201-1612	1.6mm x 12mm Partially Threaded AUTO-DRIVE™ Screw

1.6mm Standard Bone Graft Fixation Screws

208-1606	1.6mm x 6mm Screw
208-1608	1.6mm x 8mm Screw
208-1610	1.6mm x 10mm Screw
208-1612	1.6mm x 12mm Screw
208-1614	1.6mm x 14mm Screw



3885 Arapaho Road

Addison, Texas 75001

Customer Service: 800.456.7779

Main: 972.677.4600

Fax: 800.390.2620

www.osteomedcorp.com



A COLSON ASSOCIATE

*Indicated for tacking on the dental alveolar applications only.

P/N 030-1063 Rev. G