

EXTRACTOR™

Bone Marrow Aspiration Needle



ACCESSORIES

EXTRACTOR

The Extractor High Performance Bone Marrow Aspirator is engineered to be the surgeon's first choice for procedures requiring bone marrow harvesting. The longer, six-ported cannula is designed for more flexible positioning and enhanced marrow extraction. The ergonomically designed "Twist Lock" handle combines with a "Double Diamond" stylet point for greater ease of use and enhanced patient comfort.

Ideal for all Patients

- Designed for all body types
- 6 inch ported cannula keeps handle above instrumentation
- Designed for improved patient comfort

Innovative Design

- Large side ports enhance bone marrow harvesting
- The "Twin Peaks" cannula tip provides superior cutting ability
- Large 11 gauge cannula

Easy to Use

- Ergonomically designed handle allows the clinician to apply consistent pressure for greater comfort without excessive hand fatigue
- A luer lock connection enables easy and secure syringe attachment

Disposable

- Single-use patient device
- Conveniently packaged in a sterile pouch
- Provided with a protective sleeve for safe disposal



BIOSTRUCTURES LLC | T 949.553.1717 | F 949.553.0407 | 1201 Dove Street, Suite 470 • Newport Beach • CA 92660 | www.biostructures.net

