

# osteosPAN

## MORPHEUS

BONE GRAFT SUBSTITUTE

- Moldable
- Ready to use out of the package
- Osteoconductive
- Biocompatible
- Resorbable (6-12 months)
- Macroporous for bone regeneration
- Maintains graft material in operative site & resists irrigation
- Proprietary organic alkylene oxide carrier



**BIOGennix**

18011 Sky Park Circle, Suite M • Irvine, CA 92614  
info@Biogennix.com • www.Biogennix.com  
T: 949•253•0994 • F: 949•266•5800



# BIOGENNIX

Developing Unique Products for Bone Fusion

## osteoSPAN MORPHEUS

BONE GRAFT SUBSTITUTE

osteoSPAN Morpheus was evaluated in critically sized bone defects involving the metaphyseal and diaphyseal regions of the tibia where cortical and cancellous bone was present. The product was gently packed into the defect site without allograft or autograft additives. Performance was judged by micro-CT, decalcified (H&E staining) and non-decalcified (SEM-BSE) histology. Three-dimensional volumetric analyses of the defect area compared favorably with histological sections and showed Morpheus outperformed Actifuse ABX in both percent bone fraction and implant resorption.



CAT. NO.	VOL.
003-MOR	3 cc
006-MOR	6 cc
012-MOR	12 cc



Also available osteoSPAN  
blocks and granules.

# BIOGENNIX