



100% Demineralized Cancellous sponge blocks and fillers for use with intervertebral cages, arthrodesis, or in fracture sites.

Elastic and Compressible

Once hydrated, Vegagraft can compress to less than one third of its original size to precisely fill voids, intervertebral cages and irregularly shaped defects. When compressed and packed into a void, it offers maximum surface contact with the surrounding autologous structures and does not migrate from the surgical site.

Available in a vast array of sizes.

Osteoinductive

Demineralization revolutionizes bone handling properties while maintaining natural growth factors. The available BMPs in the scaffold provide osteoinductive potential while the 3 dimensional lattice's open structure and large surface area serve as an optimal natural scaffold for osteoblast proliferation and new bone growth

Ordering information:

ITEM NUMBER	PRODUCT DESCRIPTION
VX-2826	26mm x 19mm x 07 mm
VX-2830	30mm x 20mm x 05 mm
VX-2850	50mm x 20mm x 05 mm
VX-2857	50mm x 20mm x 07 mm
VX-3116	12mm x 12mm x 12mm
VX-2814	14mm x 14mm x 14mm
VX-3144	20mm x 20mm x 10mm
VX-3115	22mm x 15mm x 05mm
VX-3153	22mm x 22mm x 10mm



To place an order, please contact: customerservice@vivex.com
Vivex Biomedical, Inc., 1755 West Oak Pkwy., Suite 200, Marietta, Georgia 30062 855.MY-VIVEX www.vivex.com

MKG-30003
Rev. 00