

Amendia

The Elements of Healing

DIAMOND

ANTERIOR CERVICAL PLATE

Z_A

Z_C

Z_L

Z_O

P_x

L

S_p

S_a

Z_P

C_m

Amendia

4.0v
4.0f
4.3v
4.3f

D_d

DIAMOND
ANTERIOR CERVICAL PLATE



S_v_T

Z_T

Z_v

FEATURES AND BENEFITS

Amendia

4.0v
4.0f
4.3v
4.3f

Dd

DIAMOND
ANTERIOR CERVICAL PLATE


Diamond Plate for Anterior Cervical Fixation

Features & Benefits

- Constructed of implant grade Titanium Alloy (Ti 6Al-4V ELI)
- Plate thickness: 2.4 mm
- Maximum plate width: 19 mm
- Plate lengths: 22 mm through 92 mm (measured from actual ends of plate)
- Single level up to four level plates
- Patented Nitinol (NiTi) self-locking screw retention system
- Ø4.0 mm standard and Ø4.3 mm revision screws
- Fixed and variable screws for rigid, dynamic or hybrid stabilization
- Fixed screws allow for 1° of freedom (+/- 0.5°)
- Variable screws allow for 30° of freedom (+/- 15°)
- Screws are self-tapping and offered with blunt tips or sharp self-drilling tips
- Color coded screws based on length of screws and fixed/variable head identification
- Pre-lordosed with medial/lateral curve
- Diamond window allows for graft visualization

Available Screw Sizes

Diameter (mm)	Lengths (mm)
4.0	12,14,16,18
4.3	12,14,16,18



Amendia

The Elements of Healing

1755 West Oak Parkway

Marietta, GA 30062

Toll Free: 877-755-3329

www.amendia.com

info@amendia.com