

# Amendia

The Elements of Healing

## FORAMINAL SCOPES SYSTEM

Z<sub>C</sub>

Z<sub>A</sub>

Z<sub>L</sub>

Z<sub>O</sub>

Px

Z<sub>P</sub>

Sp

Sa

Dd

Cm

Amendia

A<sup>2</sup>

ASAP  
SPINE ENDOSCOPE SYSTEMS

Sv<sub>T</sub>

Z<sub>T</sub>

Z<sub>V</sub>

## FEATURES AND BENEFITS




# Foraminal Scopes System


## Features & Benefits

- Unique Side Cutting Manual Burrs - Safe to use for Bone Prep
- Plain and Cannulated - 3.5 mm, 5.0 mm, 6.3 mm, 6.8 mm, and 9.8 mm
- Modular Kerrisons Sizes - 3.0 mm, 3.5 mm, 4.0 mm with 90 and 135 Degrees Bites
- Cannulated Trephines
- Palpation and Irrigation Probe
- Robust Micro Spoon and Alligator Bite Graspers - 2.7 mm, 3.0 mm, 3.5 mm, and 4.0 mm
- Punches 2.7 mm and 3.5 mm

Part Number	Description
10-1496-00	Oval Foraminal Scope



Part Number	Description
10-1492-00	Round Foraminal Scope



- 6.3 mm x 6.9 mm for Ease of Access to Foramen
- 4.1 mm Working Channel for Faster Work and Decompressions
- Allows Robust 4.0 mm Instrumentation
- Dual Irrigation - One is Static, One Rotates 360° for Balance
- 1.8 mm Irrigation Channel Clearer Video Picture
- Suction on the Outer Sheath for Better Fluid Flow
- Option of 20 or 30 Degree Viewing Angles
- Pistol Grip Light Connector and Camera Head Connection
- 8.0 mm Outer Sheaths Beveled and Fenestrated

- 6.3 mm for Smaller Anatomy
- 3.7 mm Working Channel for Robust 3.5 mm Instrument Use
- Traditional Suction and Irrigation on Scope Barrel
- Options of 20 and 30 Degree Viewing Angle
- Pistol Grip Light and Camera Head Connector
- 7.0 mm Outer Sheaths - Beveled and Fenestrated

Manufactured For:



1755 West Oak Parkway  
 Marietta, GA 30062  
 Toll Free: 877-755-3329  
[www.amendia.com](http://www.amendia.com)  
[info@amendia.com](mailto:info@amendia.com)