

Amendia

The Elements of Healing

ZEUS-T

TRANSFORAMINAL LUMBAR INTERBODY

Z_A

Z_C

Z_L

Z_O

P_x

L

S_p

S_a

D_d

C_m

Amendia

Z_T

ZEUS TRANSFORAMINAL
LUMBAR INTERBODY

Sv_T

Z_P

Z_V




FEATURES AND BENEFITS




Zeus- Transforaminal Lumbar Interbody

Features & Benefits

- Constructed of PEEK
- Tantalum radiopaque markers to optimize visibility and placement
- Angled teeth prevent implant migration
- Large graft cavity to facilitate fusion
- Design ensures a high degree of compressive strength and dimensional stability
- Bullet Nose option has tapered leading edge that distracts endplates and facilitates delivery



Available Sizes (mm)	
10 W x 28 L x 7 H	11 W x 36 L x 7 H
10 W x 28 L x 8 H	11 W x 36 L x 8 H
10 W x 28 L x 9 H	11 W x 36 L x 9 H
10 W x 28 L x 10 H	11 W x 36 L x 10 H
10 W x 28 L x 11 H	11 W x 36 L x 11 H
10 W x 28 L x 12 H	11 W x 36 L x 12 H
10 W x 28 L x 13 H	11 W x 36 L x 13 H
10 W x 28 L x 14 H	11 W x 36 L x 14 H
10 W x 28 L x 15 H	11 W x 36 L x 15 H



Available Sizes (mm)
11 W x 28 L x 7 H
11 W x 28 L x 8 H
11 W x 28 L x 9 H
11 W x 28 L x 10 H
11 W x 28 L x 11 H
11 W x 28 L x 12 H
11 W x 28 L x 13 H
11 W x 28 L x 14 H
11 W x 28 L x 15 H
11 W x 36 L x 7 H
11 W x 36 L x 8 H
11 W x 36 L x 9 H
11 W x 36 L x 10 H
11 W x 36 L x 11 H
11 W x 36 L x 12 H
11 W x 36 L x 13 H
11 W x 36 L x 14 H
11 W x 36 L x 15 H