

Amendia

The Elements of Healing

SAMSON

VERTEBRAL BODY REPLACEMENT SYSTEM

Z_C

Z_A

Z_L

Z_O

Z_P

L

Sp

Px

Dd

Cm

Amendia

Sa

SAMSON CORPECTOMY
VERTEBRAL BODY REPLACEMENT



Sv_T

Z_T

Z_V

FEATURES AND BENEFITS

Samson Corpectomy Vertebral Body Replacement System

Features & Benefits

- Constructed of PEEK
- Tantalum radiopaque markers to optimize visibility and placement
- Pyramidal teeth provide multi- directional resistance to migration
- Hollow core allows for large amounts of graft material to facilitate fusion
- Cylindrical shape
- Design ensures a high degree of compressive strength and dimensional stability



14 mm Tall 32 mm Tall 48 mm Tall

Available Sizes (mm)

14 D x 14 H	14 D x 31 H
14 D x 15 H	14 D x 32 H
14 D x 16 H	14 D x 33 H
14 D x 17 H	14 D x 34 H
14 D x 18 H	14 D x 35 H
14 D x 19 H	14 D x 36 H
14 D x 20 H	14 D x 37 H
14 D x 21 H	14 D x 38 H
14 D x 22 H	14 D x 39 H
14 D x 23 H	14 D x 40 H
14 D x 24 H	14 D x 41 H
14 D x 25 H	14 D x 42 H
14 D x 26 H	14 D x 43 H
14 D x 27 H	14 D x 44 H
14 D x 28 H	14 D x 45 H
14 D x 29 H	14 D x 46 H
14 D x 30 H	14 D x 47 H
	14 D x 48 H

Amendia

The Elements of Healing

1755 West Oak Parkway

Marietta, GA 30062

Toll Free: 877-755-3329

www.amendia.com

info@amendia.com