

Amendia

The Elements of Healing

SAVANNAH-T

HIGH TOP PEDICLE SCREW

Z_A

Z_C

Z_L

Z_O

Px

L

Sp

Sa

Dd

Cm



Amendia

5.5
6.5
7.5

Sv
H

SAVANNAH-T
HIGH TOP PEDICLE SCREW

Z_P

Z_T

Z_V

FEATURES AND BENEFITS

Amendia

5.5
6.5
7.5



Sv
H

SAVANNAH-T
HIGH TOP PEDICLE SCREW

Savannah-T High Top Pedicle Screw

Features & Benefits

- Constructed of implant grade Titanium Alloy (Ti 6Al-4V ELI)
- Cannulated screws
- Self-tapping "Grip Quick" Thread®
- 5.5 mm bulleted rods - straight or pre-lordosed
- Rods offered in 5 mm increments from 30 mm through 130 mm
- Rods offered in various increments from 140 mm through 300 mm (See catalog for available sizes)
- Set screw designed to prevent tulip splaying
- 13 mm extended thread for rod reduction
- Easy release break away towers
- Minimally invasive delivery reduces soft tissue disruption
- Anti-cross thread feature in tulip and set screw
- Connected Tabs or Open Top options available



Available Sizes	
Diameter (mm)	Lengths (mm)
5.5	30-50
6.5	30-65
7.5	30-65

Amendia

The Elements of Healing

1755 West Oak Parkway

Marietta, GA 30062

Toll Free: 877-755-3329

www.amendia.com

info@amendia.com