

# Amendia

The Elements of Healing

## SPARTAN S<sup>3</sup>

FACET SCREW SYSTEM

Z<sub>A</sub>

Z<sub>C</sub>

Z<sub>L</sub>

Z<sub>O</sub>

Px

L

Z<sub>P</sub>

Sa

Dd

Cm

Amendia

4.3  
5.0

Sp

SPARTAN  
FACET SCREW

Sv<sub>T</sub>

Z<sub>T</sub>

Z<sub>V</sub>

## FEATURES AND BENEFITS

## Spartan Facet Screw

## Features &amp; Benefits

- Constructed of implant grade Titanium Alloy (Ti 6Al-4V ELI)
- “Grip Quick” Thread® allows for faster screw delivery
- Spiral self-tapping flutes reduce drive torque and improve integration
- Lag screws have increased shank diameter for stability and strength
- Anti back-out threads under flared head provide rapid change in torque resistance when seated
- Optional skirt provides force distribution to minimize risk of damage to the facet
- Robust 1.4 mm guide wire resists kinking
- Threaded guide wire reduces effort, is easy to control without rapid plunge and anchors into the facet
- Drill marked for depth relation
- Driver retains screw for retraction or revision
- Internal square drive cannot disengage with retainer collar in place
- Ratcheting silicone axial driver handle
- Driver disassembles for thorough decontamination and sterilization

Diameter (mm)	Length (mm)	Skirt Option
4.3	20	No
4.3	25	No
4.3	30	No
4.3	35	No
5.0	25	Yes
5.0	30	Yes
5.0	35	Yes
5.0	40	Yes



**Amendia**  
The Elements of Healing

1755 West Oak Parkway  
Marietta, GA 30062  
Toll Free: 877-755-3329  
[www.amendia.com](http://www.amendia.com)  
[info@amendia.com](mailto:info@amendia.com)