

Proprietary Innovations

Medline revolutionized latex-free surgical gloves with SensiCare, the first polyisoprene glove on the market. Today, Medline's SensiCare gloves are being used by some of the most prestigious healthcare facilities worldwide. In addition, Medline is the only surgical glove supplier that manufactures the raw material polyisoprene, allowing complete quality control from start to finish.

Choose SensiCare and Never Go Back

Whether your facility is ready to go completely latex-free or just interested in upgrading to an exceptional synthetic material, Medline has a comprehensive line of SensiCare polyisoprene gloves with advanced technologies to satisfy even the most demanding surgeons.

Experience a "BE FREE DAY"

To create awareness about the risks associated with latex, Medline offers a "BE FREE DAY". At NO cost to your facility Medline will provide a day's worth of our SensiCare surgical gloves as a no-pressure opportunity to introduce your staff to the newest latex-free technology. Contact your Medline sales representative or go to medline.com/gloves for more information.



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)
www.medline.com
info@medline.com

Medline Canada
1-800-396-6996
www.medline.ca
canada@medline.com

Medline México
01-800-831-0898
www.medlinemexico.com
mexico@medline.com

Some products may not be available for sale in Mexico or Canada.
We reserve the right to correct any errors that may occur within this sell sheet.
©2011 Medline Industries Inc. All rights reserved. Isolex is a trademark and
Medline and SensiCare are registered trademarks of Medline Industries, Inc.

MKT1330977 / LIT102R / 5M / DLP4



SENSICARE® LATEX-FREE SURGICAL GLOVES
Advanced Performance and Protection



Isolex™: SensiCare's Breakthrough Technology

Medline's SensiCare surgical gloves are made from Isolex, a proprietary synthetic polyisoprene. This material has a molecular structure virtually identical to natural rubber latex but without the harmful latex proteins. As a result, SensiCare surgical gloves are softer, more elastic, and more comfortable than latex to satisfy clinical needs and support safety initiatives.



SensiCare Glove	Ideal Use	Details	Item Numbers
SensiCare® with Aloe Vera	All general surgery procedures	The First Polyisoprene Glove 10 mil thickness, smooth, white, half sizes 5.5–9.0	MSG1055 to MSG1090, 25/bx, 100/cs
SensiCare® LT with Aloe Vera	All general surgery procedures	Enhanced Grip, Superior Dexterity 10 mil thickness, textured, white, half sizes 5.5–9.5*	MSG1155 to MSG1195, 25/bx, 100/cs
SensiCare® SLT	All general surgery procedures	Improved Fit and Sensitivity 9.5 mil thickness, textured, cream, half sizes 5.5–9.0	MSG1555 to MSG1590, 50/bx, 200/cs
SensiCare® Green with Aloe Vera	Double-gloving in cardiac, trauma, and high-risk procedures	Microthin for Easy Double-Gloving 8.5 mil thickness, smooth, green, half sizes 5.5–9.0	MSG1255 to MSG1290, 25/bx, 100/cs
SensiCare® Micro	Neuro, vascular, and cardiac procedures	Offering a High Level of Tactility 8.5 mil thickness, smooth, white, half sizes 5.5–9.0	MSG1655 to MSG1690, 50/bx, 200/cs
SensiCare® Ortho	Orthopaedic procedures	Thicker for Added Protection 14 mil thickness, textured, white, half sizes 6.0–9.0	MSG1460 to MSG1490, 25/bx, 100/cs
SensiCare® LT Custom Fit	All general surgery procedures	Tailored for Hard-to-Fit Hands 10 mil thickness, smooth, white, half sizes 7.0–8.5	MSG1170C to MSG1185C, 25/bx, 100/cs
Sensicare® Shield	Fluoroscopic procedures to shield hand from scattered radiation	Radiation Attenuation 11.8 mil thickness, textured, black, half sizes 6.5–9.0	MSG3965 to MSG3990, 5/bx, 20/cs

*Medline is the only major surgical glove company in the U.S. offering size 9.5