







Apollo Medial Suture Anchor Delivery System

Features and Benefits



Apollo Medial Suture Anchor Delivery System

Features and Benefits

- Preloaded fully disposable delivery system with fully contained sutures
- Two piece 'Stable Drive' handle design
- PEEK anchor available in 4.5mm, 5.5mm, and 6.5mm diameters
- Integrated stainless steel punch for insertion without predrilling
 - Stainless steel punch removed upon insertion
- Double lead tapered threads for faster delivery
- Two #2 High Strength UHMWPE/Polyester sutures
 - One White suture
 - One White with Green accent suture

