






Apollo Lateral Anchor Delivery System

Features and Benefits



Apollo Lateral Anchor Delivery System

Features and Benefits

- Preloaded fully disposable delivery system with integrated stainless steel punch
- Two piece 'Stable Drive' handle design
- PEEK anchor available in 4.5mm, 5.5mm, and 6.5mm diameters
- Dual slotted Suture Bullet for suture capture from multiple locations and knotless repairs
 - Captures up to eight #2 sutures
 - Delivered to 20mm depth
- Integrated stainless steel punch for insertion with or without predrilling
 - Punch locked for suture bullet delivery and graft tensioning
 - Punch unlocked for anchor delivery
 - Removed upon full delivery
- Double lead tapered threads for faster delivery and superior pullout strength
- Designed for maximum suture and anchor engagement for superior resistance to pull out

