

## PinkCare PF Nitrile MG505 Series



## Characteristics and Overview:

Medgluv latex-free exam glove, PinkCare Nitrile will offer you the protection of Nitrile and feeling of latex. Featured with an advanced stretch formulation, it will provide you ultimate comfort for extended periods. These latex-free gloves offer a high valued and dependable solution to those allergic to natural latex protein. These gloves, formulated with 100% synthetic nitrile polymer, are latex free, powder free and odor free. The unique thermal sensitive polymer compound absorbs body heat and easily conforms to user's hands; thus, everts virtually no pressure on the hands and reduces fatigue after extended wear.

### Sizes

MG555XS	MG555S	MG555M	MG555L	MG555XL
X-Small	Small	Medium	Large	X-Large
200 ea/bx	200 ea/bx	200 ea/bx	200 ea/bx	200 ea/bx



## Glove Specifications

	ASTM	PinkCare
Length (mm)	Min 220	240
Palm width (mm)	95 ± 10	95 ± 10
Cuff Thickness (mm)	N/A	0.06 (2.4 mils)
Palm Thickness (mm)	min. 0.05	0.07 (2.8 mils)
Finger Thickness (mm)	min. 0.05	0.09 (3.5 mil)
Tensile Strength (MPa)	min. 14	29
Ultimate Elongation (%)	min 400	570

## Glove Attributes

Material	Nitrile
Sterility	Non-Sterile
Color	Pink
Sizes	X-Small, Small, Medium, Large, X-Large
Glove Interior	Powder-Free, Chlorinated
Glove Exterior	Textured Fingers- tips
Packaging	200 gloves per box, 10 boxes per case

## Quality Standards

Exceeds current ASTM 6319 standard for critical defects (AQL 2.5). AQL for critical defects is 1.5.  
Manufactured in accordance with Quality System ISO 9001.

**\*\* Tested for use with chemotherapy drugs using ASTM Standard D6978**

**\*\* PinkCare nitrile gloves symbolizes Medgluv's commitment to supporting breast cancer research and outreach. Medgluv will donate 5% of PinkCare Nitrile Gloves profits toward breast cancer research and awareness programs**