

NitriSkin Sterile MGSE502 & MGSE552 Series



Characteristics and Overview:

NitriSkin Sterile Singles and Pairs are now the ultimate choice use of medical professionals in high demand environments like Obstetricians and Gynecologists, intensive care unit ICC/CCU, traumas, burn centers and wound care that requires a high quality latex-free sterile glove that delivers maximum barrier protection.

Sizes

MGSE502S Small 100 ea/bx	MGSE502M Medium 100 ea/bx	MGSE502L Large 100 ea/bx	MGSE502XL X-Large 100 ea/bx
MGSE552S Small 50 pr/bx	MGSE552M Medium 50 pr/bx	MGSE552L Large 50 pr/bx	MGSE552XL X-Large 50 pr/bx

Glove Specifications

	ASTM	NitriSkin Sterile
Length (mm)	Min 220	240
Palm width (mm)	95 ± 10	95 ± 10
Cuff Thickness (mm)	N/A	0.12 (5.0 mil)
Palm Thickness (mm)	min. 0.09	0.14 (5.5 mil)
Finger Thickness (mm)	min. 0.05	0.15 (5.8 mil)
Tensile Strength (MPa)	min. 14	25
Ultimate Elongation (%)	min 400	600

Glove Attributes

Material	Nitrile
Sterility	Sterile
Color	Blue
Sizes	Small, Medium, Large, X-Large
Glove Interior	Powder-Free, Chlorinated
Glove Exterior	Textured
Packaging	Singles - 100 gloves per box, 4 boxes per case Pairs - 50 pairs per box, 4 boxes per case

Quality Standards

Exceeds current ASTM D6319 standard for critical defects (AQL 2.5). AQL for critical defects is 1.5.
Manufactured in accordance with Quality System ISO 9001.

**** Tested for use with chemotherapy drugs using ASTM Standard D6978**