

NitraCare Plus PF Nitrile MNE505 Series



Characteristics and Overview:

Medgluv best selling latex-free glove, NitraCare Nitrile will offer you the protection of Nitrile and feeling of latex. Featured with an advanced stretch formulation, it will provide you ultimate comfort for extended periods. These latex-free gloves offer a high valued and dependable solution to those allergic to natural latex protein. These gloves, formulated with 100% synthetic nitrile polymer, are latex free, powder free and odor free. The unique thermal sensitive polymer compound absorbs body heat and easily conforms to user's hands; thus, everts virtually no pressure on the hands and reduces fatigue after extended wear.

Sizes

MNE5051	MNE5052	MNE5053	MNE054	MNE5054	MNE5055
X-Small	Small	Medium	Large	X-Large	2X-Large
100 ea/bx	100 ea/bx	100 ea/bx	100 ea/bx	100 ea/bx	90 ea/bx

Glove Specifications

	ASTM	NitraCare
Length (mm)	Min 220	240
Palm width (mm)	95 ± 10	95 ± 10
Cuff Thickness (mm)	N/A	0.08 (2.8 mil)
Palm Thickness (mm)	min. 0.05	0.08 (3.1 mil)
Finger Thickness (mm)	min. 0.05	0.09 (3.6 mil)
Tensile Strength (MPa)	min. 14	29
Ultimate Elongation (%)	min 400	570

Glove Attributes

Material	Nitrile
Sterility	Non-Sterile
Color	Cobalt Blue
Sizes	X-Small, Small, Medium, Large, X-Large, 2X-Large
Glove Interior	Powder-Free, Chlorinated
Glove Exterior	Textured Fingers- tips
Packaging	100 gloves per box, 10 boxes per case 2XL packed 90 gloves per box

Quality Standards

Exceeds current ASTM 6319 standard for critical defects (AQL 2.5). AQL for critical defects is 1.5.

Manufactured in accordance with Quality System ISO 9001

**** Tested for use with chemotherapy drugs using ASTM Standard D6978**