

## XP Gold PF Latex MG1210 Series



### Characteristics and Overview:

XP Gold Powder-Free Latex Exam glove is thicker and longer than ordinary exam gloves for extra protection in high risk situations, yet providing exceptionally good tactile sensitivity. Ultra low modulus compound and polymer coated for comfortable use for extended periods.

XP Gold glove meets or exceeds the emergency medical glove requirements of NFPA 1999, 2008 edition.

A leading High-Risk EMS glove, with a low drag textured finish ideal for Police and Corrections applications.

#### Sizes

MG1210S	MG1210M	MG1210L	MG1210XL	MG1210XXL
Small	Medium	Large	X-Large	2X-Large
50 ea/bx	50 ea/bx	50 ea/bx	50 ea/bx	44 ea/bx

### Glove Specifications

	ASTM	XP Gold
Length (mm)	Min 280	290
Palm width (mm)	95 ± 10	95 ± 10
Cuff Thickness (mm)	N/A	0.20 (7.9 mil)
Palm Thickness (mm)	min. 0.05	0.22 (8.5 mil)
Finger Thickness (mm)	min. 0.05	0.25 (10.0 mil)
Tensile Strength (MPa)	min. 14	38
Ultimate Elongation (%)	min 400	600

### Glove Attributes

Material	Latex
Sterility	Non-Sterile
Color	Beige
Sizes	Small, Medium, Large, X-Large, 2X-Large
Glove Interior	Powder-Free, Chlorinated
Glove Exterior	Full Textured
Packaging	50 gloves per box, 10 boxes per case 2XL packed 42 gloves per box



### Quality Standards

Exceeds current ASTM D3578 standard for critical defects (AQL 2.5). AQL for critical defects is 1.5.

Manufactured in accordance with Quality System ISO 9001

Quality sampled in accordance with MIL.STD 105E

**\*\* Tested for use with chemotherapy drugs using ASTM Standard D6978**

**THIS GLOVE MEETS THE EMERGENCY MEDICAL GLOVE REQUIREMENTS OF NFPA 1999, STANDARD ON PROTECTIVE CLOTHING FOR EMERGENCY MEDICAL OPERATIONS, 2008 EDITION.**

