



CATHY® IV CATHETERS

A New Level of Performance

Cathy® IV Catheters provide you with catheters that are easy to use with performance attributes that reduce peelback and improve patient comfort. Cathy® IV Catheters provide the performance and quality you need while improving user convenience and patient measures.



Bringing **Reliability, Innovation and Value** to Healthcare®

Cathy® IV Catheters

The Performance and Convenience You Require

- Siliconized short bevel needle and three facet cut needle point provides smoother insertion
- Catheter manufactured from proprietary polymer
- Proprietary polymer rigid at room temperature reducing peelback and softens inside body for improved patient comfort
- Proprietary polymer is radiopaque
- Color coded for easy gauge identification
- Catheter contains angled, grooved wings which improve stabilization and provide for easy placement


The Highest Quality Levels

- Fully automated assembly and 100% vision camera inspection assures proper catheter tip geometry
- Computerized grinding and 100% vision camera inspection ensures sharp, consistent, high quality needles

The Sizes You Demand

- We offer 22G and 24G sizes for your users

Cathy® IV Catheters Needle Gauge Color Chart



To Inquire about MYCO's Products or receive complimentary evaluation samples, please contact:



1.800.454.6926 or 1.919.460.2535

Info@mycomedical.com



Bringing **Reliability, Innovation and Value** to Healthcare®