

MICROFUSE RAPID RATE SYRINGE INFUSER

THE MICROFUSE RAPID RATE SYRINGE INFUSER IS DESIGNED FOR **SPECIALTY APPLICATIONS INFUSING DOSES QUICKLY.**



64004FU-0010 Rev C



SPECIFICATIONS

Size: 9.6" x 2.7" x 4"

Weight: 25 oz. with batteries

Compatible Syringes: BD, Monoject and Terumo

Syringe Sizes: 30 – 60 cc

Accuracy of Flow Rate: +/- 3%*

Power: 2, "C" cell alkaline batteries

Battery Life: Typically 400 operational hours

Infusion Rates: 1 to 6 minutes for full syringe
2 to 8.4 minutes for full syringe

Alarms: Occlusion, end of infusion, low battery and internal malfunction

* excluding syringe variations



an infusion of new ideas

NUMIA
MEDICAL TECHNOLOGY, LLC.

Numia Medical Technology, LLC | 84 Farrant Street | Newport, Vermont 05855
tel 802 323-0101 | fax 802 334-6971 | info@numiamedical.com | www.numiamedical.com