

MICROFUSE DUAL RATE SYRINGE INFUSER

THE MICROFUSE DUAL RATE SYRINGE INFUSER IS DESIGNED FOR INFUSION OF **NON-RATE CRITICAL INTRAVENOUS DRUGS**.



6400-IFU-0008 Rev E



SPECIFICATIONS

Size: 9.6" x 2.7" x 4"

Weight: 25 oz. with batteries

Compatible Syringes: BD, Monoject and Terumo

Syringe Sizes: 5 – 140mL

Accuracy of Flow Rate: +/- 3%*

Power: 2, "C" cell alkaline batteries

Battery Life: Typically 400 operational hours

Infusion Rates: Normal: 20 – 40 minutes, Slow: 60 – 120 minutes

Alarms: Occlusion, end of infusion, low battery and internal malfunction

* excluding syringe variations



an infusion of new ideas

NUMIA
MEDICAL TECHNOLOGY, LLC.

Numia Medical Technology, LLC | 84 Farrant Street | Newport, Vermont 05855
tel 802 323-0101 | fax 802 334-6971 | info@numiamedical.com | www.numiamedical.com