



BreathCare BiPAP Device  
**Favor a Better Life**

# Sleep Fairy

- Advanced high precision 3D gyroscope technology.
- 8 parameters: SpO2, Pulse, Airflow, Snoring, Chest & Abdomen Movement, Body position, Body movement, CPAP Pressure.
- 2.8 inch TFT color display.
- Easy operation.
- The battery can support over 20hours data recording.
- Large capacity rechargeable lithium battery.
- Patented chest & abdomen signal acquisition.
- Intelligent analysis software.
- High speed SD card for data record without computer operation limited.





## **Comfort:**

EPR(Expiration Pressure Relief) ensure comfort degree  
Automatic Synchronization Technology ensure the ventilator  
and patient to achieve perfect synchronization,  
which make patients nice experience

## **Quiet : Noise $\leq 28\text{dB}$**

High efficiency brushless DC motor, 10MS of rapid response  
Patent for cancel muffling technology

## **Easy Operation**

Special water tank design to make operation very easy and  
humanized





One or two cups water per night, anti counter-current, Anti-dry safety protection technology



Humidifier can heat directly without tank



Up and Down Key makes operation very easy and durable



Real-time display treatment data

## Auto CPAP (Silver)

**Working Mode:** CPAP, Auto

**Sound:** <28dB

**Screen:** 2.8 inch TFT, 320 x 240

**Pressure Range:** 4-20cmH<sub>2</sub>O, auto adjust

**Ramp Item:** 0-60min adjustable 1min/Step

**Water Capacity:** 220ML max

**Humidifier Level:** 0-5 level

**Humidifier:** Integrated Heated Humidifier

**Date Storage:** SD card

Auto CPAP mode can work for both OSA and CSA

Leak Compensation

Auto-Altitude Compensation

## BiPAP (Grey)

**Working Mode:** CPAP, Auto CPAP, S, T, ST

**Sound:** <28dB

**Screen:** 2.8 inch TFT, 320 x 240

**Pressure Range:** 4-25cmH<sub>2</sub>O(1);  
4-30cmH<sub>2</sub>O(2)

**IPAP:** 4-25cmH<sub>2</sub>O (1); 4-30cmH<sub>2</sub>O(2)

**EPAP:** 4-20 cmH<sub>2</sub>O

**Ramp Item:** 0-60min adjustable 1min/Step

**Humidifier:** Integrated Heated Humidifier

**Date Storage:** SD card

Leak Compensation

Auto-Altitude Compensation