

**Material Safety Data Sheet**  
According to Regulation (EC) No: 1907/2006

**Dispette**  
Revision Date: June 2013

**1. IDENTIFICATION OF THE SUBSTANCES/PREPARATION AND OF THE COMPANY UNDERTAKING**

PRODUCT NAME: Dispette

APPLICATION: Disposable tube for erythrocyte sedimentation rate estimation

MANUFACTURER: Guest Scientific A G  
Gewerbestrasse 10  
CH-6330 Cham  
Switzerland

☎ ++41(0)41 741 4748  
Fax ++41(0)41 741 3966

email: guest.scientific@bluewin.ch

**2. HAZARDS IDENTIFICATION**

No hazards identified - no chemical constituents

**3. COMPOSITION/INFORMATION ON INGREDIENTS**

Pipette manufactured from clear polystyrene. Printed with graduations 200 to 0 mm. Filling reservoir manufactured from blue polypropylene.

**4. FIRST AID MEASURES**

Not applicable

**5. FIRE-FIGHTING MEASURES**

No special measures required

**6. ACCIDENTAL RELEASE MEASURES**

Not applicable

**7. HANDLING AND STORAGE**

Store in a cool dry place

**8. EXPOSURE CONTROLS/PERSONAL PROTECTION**

WEL = Workplace Exposure Limit. Not applicable

**9. PHYSICAL AND CHEMICAL PROPERTIES**

APPEARANCE clear polystyrene tube with white printed graduations & blue filling reservoir

ODOUR Not applicable

SOLUBILITY Insoluble in water

pH-VALUE, DILUTED SOLUTION Not applicable

DECOMPOSITION TEMPERATURE Polystyrene 240°C. Polypropylene between 130°C & 171°C

#### 10. STABILITY AND REACTIVITY

STABILITY: Stable under normal temperature conditions.

CONDITIONS TO AVOID: Avoid contact with solvents. Avoid exposure to high temperatures

HAZARD DECOMPOSITION PRODUCTS: Complete oxidation of polystyrene produces carbon dioxide and water, Partial oxidation may release carbon monoxide.

#### 11. TOXICOLOGICAL INFORMATION

Not applicable

#### 12. ECOLOGICAL INFORMATION

DEGRADABILITY: This product is not biodegradable.

#### 13. DISPOSAL CONSIDERATIONS

DISPOSAL METHODS: Dispose of waste and residues in accordance with local authority requirements.

#### 14. TRANSPORT INFORMATION

UK ROAD CLASS NA

PROPER SHIPPING NAME NA

ROAD TRANSPORT NOTES NA

AIR TRANSPORT NOTES Not Classified

#### 15. REGULATORY INFORMATION

Polystyrene is classified according to DIN4102 as a B3 product – meaning highly flammable or easily ignited.

No Special labelling requirements

No Risk phrases

No Safety Phrases

#### 16. OTHER INFORMATION

##### GENERAL INFORMATION

REVISION DATE 14 06 2013

Rev. No./REPL. SDS GENERATED 4

#### DISCLAIMER

This information relates only to the specific materials designated and may not be valid for such material used in combination with any other materials or in any process. Such information is, to the best of the company's knowledge and belief, accurate and reliable as of the date indicated. However, no warranty, guarantee or representation is made to its accuracy, reliability or completeness. It is the user's responsibility to satisfy himself as to the suitability of such information for his own particular use.

MSDS/ESR/1/Disp rev.06.13