



Comforter™ Splint



The softly padded Comforter™ Splint provides the rest needed to reduce inflammation and morning stiffness cause by Rheumatoid Arthritis.

This comfortable night splint is also used for stroke patients and those with ulnar deviation of the wrist and MP joints.

Features

- The Comforter Splint holds the wrist and MP joints of the fingers straight to reduce stress on muscles and tendons, making it a great night splint for arthritis sufferers.
- Flexible padded stays protect the fingers to lessen the damaging effects of arthritis and certain neurological impairments.
- The small joints of the fingers and the thumb are free to allow the wearer to manage bedding and nighttime activities. The wrist and the metacarpal joints are supported in a neutral position.
- The neoprene straps include finger holes for independent application and removal of the splint even by those who have difficulties with fine motor skills.
- Quilted cover is removable and machine washable. Latex-free.

Order Information

Measure width from the 2nd to the 5th MP (large knuckles) to determine size.

	Right	Left	Fits
Sm	P3002-R2	P3002-L2	2 3/4" to 3"
Med	P3002-R3	P3002-L3	3" to 3 3/8"
Lrg	P3002-R4	P3002-L4	3 3/8" to 3 5/8"



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com