



# ThumSaver™ CMC Long



The ThumSaver™ CMC Long is a full length thumb splint that protects the wrist and thumb from the pain of severe tendinitis or CMC arthritis or following scaphoid fractures while allowing for light activity.

## Features

- Carefully designed to limit thumb and wrist motion for relief from the pain of arthritis or tendinitis or for protection following cast removal.
- The ThumSaver CMC Long holds the thumb in a functional position that allows the wearer to grasp objects to perform light daily activities.
- Thin, lightweight splint includes felt pads to provide cushioning where needed.
- Seamless, molded design will not break or change shape with heat and resists stains.

## Order Information

Measure the length from the IP joint (knuckle at the end of the thumb) to the wrist. Measure around the IP joint.

	Right	Left	Length	Circumference
Sm	P3000-R2	P3000-L2	3" to 3 1/2"	2 1/2" to 2 3/4"
Med	P3000-R3	P3000-L3	3 1/2" to 4"	2 3/4" to 3 1/8"
Lrg	P3000-R4	P3000-L4	4" to 4 1/2"	3 1/8" to 3 1/4"



**Call 888-378-7763**

Mon-Fri 8:30am - 5:00pm EST for ordering details

**[www.3pointproducts.com](http://www.3pointproducts.com)**