



3pp[®] ThumSling[®] Long



The 3pp[®] ThumSling[®] Long is a lightweight, foam-lined splint that provides light support for both the wrist and thumb to reduce pain caused by thumb arthritis, ligament sprains or tendinitis. The ThumSling Long provides adjustable compression and support without restricting motion.

Features

- Supports the CMC joint (base joint) of the thumb and wrist with an adjustable wrap-around strap. Relieves pain and discomfort from thumb arthritis, tendinitis or overuse injuries.
- Recommended for those with pain primarily during activity, the 3pp ThumSling Long is not bulky or hot and is easy-to-wear at work and during leisure activities.
- Soft foam lining is comfortable against the skin and helps keep the splint in place.
- The flexible stay (included) can be added wherever extra support is needed.
- Hand or machine washable. Latex-free.

Order Information

Measure around the large knuckles of the hand. One size fits most.

	Right	Left
One Size	P3026-R	P3026-L



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com