



3pp[®] Ultra Spica



The 3pp[®] Ultra Spica is semi-flexible, weighs just one ounce and provides light, effective thumb and wrist support. Use the 3pp Ultra Spica to reduce the pain of arthritis, tendinitis or conditions caused by repetition and overuse.

Features

- Lightweight plastic is lined in soft felt for cushioning and comfort and the flexible material contours to your hand as it is strapped on.
- Allows wearer to carry out daily activities comfortably and may be used at night to rest the wrist and thumb.
- Supports the thumb and wrist in a neutral position to rest inflamed ligaments and tendons.

Order Information

Measure the length from the IP joint (knuckle at the end of the thumb) to the wrist.

	Right	Left	Fits
Sm	P3100-R2	P3100-L2	3" to 3 1/2"
Med	P3100-R3	P3100-L3	3 1/2" to 4"
Lrg	P3100-R4	P3100-L4	4" to 4 1/2"



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com