

## TECHNICAL DATA SHEET

### Conforming Bandage Roll Items # A205 to D606

#### Intended Use

Conforming bandages used to cover and protect wounds from further injuries and contamination. This type of product fits all parts of the body without restrictive movement. To be used in health care facilities.

Characteristic	Specification			
Description	Roll of loosely weaved, sheer, self-finished edge (selvage) fabric			
Material	100% Polyester. Latex-Free			
Color	White			
Dimension	Size	Item Code	Width (cm)	Length (m)
	2" x 4,1 yd	A205, A206, D205, D206	5,0	3,7
	3" x 4,1 yd	A305, A306, D305, D306	7,5	
	4" x 4,1 yd	A405, A406, D405, D406	10,0	
	6" x 4,1 yd	A605, A606, D605, D606	15,0	
Packaging	Item Code	Bandages/Box	Boxes/Case	Bandages/Case
	A205, A305, A405, D205, D305, D405	12	8	96
	A605, D605	6		48
	Item Code	Roll/Pouch	Pouches/Bag	Bags/Case
	A206, D206, A306, D306, A406, D406,	1	12	8
	A606, D606			4

Item Code	Description	Color	Specific Attributes
A205	AMD-Ritmed Conforming Stretch Bandage, 2" x 4,1 yds (5,1 cm x 3,7 m), 12 Bandages/Box, 8 Boxes/Case, 96 Bandages/Case	White	Non-Sterile
A206	AMD-Ritmed Conforming Stretch Gauze Bandage, 2" x 4,1 YD (5,1 cm x 3,7 m), 1 Roll/Pouch, 12 Pouches/Bag, 8 Bags/Case, 96 Rolls/Case	White	Individually Wrapped, Sterile
A305	AMD-Ritmed Conforming Stretch Bandage, 3" x 4,1 yds (7,6 cm x 3,7m), 12 Bandages/Box, 8 Boxes/Case, 96 Bandages/Case	White	Non-Sterile
A306	AMD-Ritmed Conforming Stretch Gauze Bandage, 3" x 4,1 YD (7,6 cm x 3,7 m), 1 Roll/Pouch, 12 Pouches/Bag, 8 Bags/Case, 96 Rolls/Case	White	Individually Wrapped, Sterile
A405	AMD-Ritmed Conforming Stretch Bandage, 4" x 4,1 yds (10,2 cm x 3,7 m), 12 Bandages/Box, 8 Boxes/Case, 96 Bandages/Case	White	Non-Sterile
A406	AMD-Ritmed Conforming Stretch Gauze Bandage, 4" x 4,1 YD (10,2 cm x 3,7 m), 1 Roll/Pouch, 12 Pouches/Bag, 8 Bags/Case, 96 Rolls/Case	White	Individually Wrapped, Sterile
A605	AMD-Ritmed Conforming Stretch Bandage, 6" x 4,1 yds (15,2 x 3,7m), 6 Bandages/Box, 8 Boxes/Case, 48 Bandages/Case	White	Non-Sterile
A606	AMD-Ritmed Conforming Stretch Gauze Bandage, 6" x 4,1 YD (15,2 cm x 3,7 m), 1 Roll/Pouch, 12 Pouches/Bag, 4 Bags/Case, 48 Rolls/Case	White	Individually Wrapped, Sterile

## TECHNICAL DATA SHEET

### Conforming Bandage Roll Items # A205 to D606

Item Code	Description	Color	Specific Attributes
D205	Distech Conforming Bandages, 2" x 4,1 yds (5,1 cm x 3,7 m), 12 Bandages/Box, 8 Boxes/Case, 96 Bandages/Case	White	Non-Sterile
D206	Distech Conforming Bandages, 2" x 4,1 YD (5,1 cm x 3,7 m), 1 Roll/Pouch, 12 Pouches/Bag, 8 Bags/Case, 96 Rolls/Case	White	Individually Wrapped, Sterile
D305	Distech Conforming Bandages, 3" x 4,1 yds (7,6 cm x 3,7m), 12 Bandages/Box, 8 Boxes/Case, 96 Bandages/Case	White	Non-Sterile
D306	Distech Conforming Bandages, 3" x 4,1 YD (7,6 cm x 3,7 m), 1 Roll/Pouch, 12 Pouches/Bag, 8 Bags/Case, 96 Rolls/Case	White	Individually Wrapped, Sterile
D405	Distech Conforming Bandages, 4" x 4,1 yds. (10,2 cm x 3,7 m), 12 Bandages/Box, 8 Boxes/Case, 96 Bandages/Case	White	Non-Sterile
D406	Distech Conforming Bandages, 4" x 4,1 YD (10,2 cm x 3,7 m), 1 Roll/Pouch, 12 Pouches/Bag, 8 Bags/Case, 96 Rolls/Case	White	Individually Wrapped, Sterile
D605	Distech Conforming Bandages – 6" x 4,1 yds (15,2 x 3,7m), 6 Bandages/Box, 8 Boxes/Case, 48 Bandages/Case	White	Non-Sterile
D606	Distech Conforming Bandages, 6" x 4,1 YD (15,2 cm x 3,7 m), 1 Roll/Pouch, 12 Pouches/Bag, 4 Bags/Case, 48 Rolls/Case	White	Individually Wrapped, Sterile

#### Regulatory

FDA Class 1, FDA Product Code MHW  
Health Canada Class 1

#### Storage Condition

Store under ambient temperature and dry conditions