

TECHNICAL DATA SHEET

Gauze Sponge Items # D-2120 to D-4800

Intended Use

May be used for one or more of the following purposes: cleansing wounds or skin, absorbing biological fluids or support organs and tissue during surgical procedures

Characteristic	Specification				
Description	Folded woven gauze sponges				
Material	100% cotton, Latex-free				
Color	White				
Ply	Item Code	Folding Method ⁽¹⁾		Ply	
	D-2800, D-3001, D-3003, D-3005, D-3800, D-4800	C-Fold		8	
	D-3002, D-3004, D-3006, D-2120, D-2212, D-3120, D-3312, D-4120, D-4412	C-Fold		12	
	D-3007, D-4416	C-Fold		16	
Dimension	Item Code	Width (cm)		Length (cm)	
	D-2800, D-2120 D-2212, D-3001 D-3002	4,6 (1,8")		5.1 (2.0")	
	D-3003, D-3004 D-3800 D-3120 D-3312	7,6 (3,0")		7,6 (3,0")	
	D-3005, D-3006, D-3007 D-4800 D-4120, D-4412 D-4416	10,2 (4.0")		10,2 (4.0")	
Packaging	Item Code	Sponges / Bag		Bags / Case	Sponges / Case
	D-3006, D3007	200		10	2000
	D-3003, D-3004, D-3005			20	4000
	D-3001			25	5000
	D-3002			40	8000
	D-2212, D-3312	Sponges / Pouch	Pouches / Box	Boxes / Case	2400
	D-4412	1	100	24	
	D-2120, D-2800	2	50	30	3000
	D-3120, D-3800		40		2400
	D-4800, D-4120		25	24	1200
	D-4416			16	800

⁽¹⁾ Refer to folding diagram


TECHNICAL DATA SHEET

Gauze Sponge Items # D-2120 to D-4800

Item Code	Description	Color	Specific Attribute
D-2120	AMD-Ritmed Gauze Sponges, 2"x2", 12 Ply, 2 Sponges/Pouch, 50 Pouches/Boxes, 30 Boxes/Case, 3000 Sponges/Case	White	Sterile
D-2212	AMD-Ritmed Gauze Sponges, 2"x2", 12 Ply, 1 Sponge/Pouch, 100 Pouches/Box, 24 Boxes/Case, 2400 Sponges/Case	White	Sterile
D-2800	AMD-Ritmed Gauze Sponges, 2"x2", 8 Ply, 2 Sponges/Pouch, 50 Pouches/Box, 30 Boxes/Case, 3000 Sponges/Case	White	Sterile
D-3001	AMD-Ritmed Gauze Sponges, 2"x2", 8 Ply, 200 Sponges/Bag, 25 Bags/Case, 5000 Sponges/Case	White	Non-Sterile
D-3002	AMD-Ritmed Gauze Sponges, 2"x2", 12 Ply, 200 Sponges/Bag, 40 Bags/Case, 8000 Sponges/Case	White	Non-Sterile
D-3003	AMD-Ritmed Gauze Sponges, 3"x3", 8 Ply, 200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	White	Non-Sterile
D-3004	AMD-Ritmed Gauze Sponges, 3"x3", 12 Ply, 200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	White	Non-Sterile
D-3005	AMD-Ritmed Gauze Sponges, 4"x4", 8 Ply, 200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	White	Non-Sterile
D-3006	AMD-Ritmed Gauze Sponges, 4"x4", 12 Ply, 200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	White	Non-Sterile
D-3007	AMD-Ritmed Gauze Sponges, 4"x4", 16 Ply, 200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	White	Non-Sterile
D-3120	AMD-Ritmed Gauze Sponges, 3"x3", 12 Ply, 2 Sponges/Pouch, 40 Pouches/Box, 30 Boxes/Case, 2400 Sponges/Case	White	Sterile
D-3312	AMD-Ritmed Gauze Sponges, 3"x3", 12 Ply, 1 Sponge/Pouch, 100 Pouches/Box, 24 Boxes/Case, 2400 Sponges/Case	White	Sterile
D-3800	AMD-Ritmed Gauze Sponges, 3"x3", 8 Ply, 2 Sponges/Pouch, 40 Pouches/Box, 30 Boxes/Case, 2400 Sponges/Case	White	Sterile
D-4120	AMD-Ritmed Gauze Sponges, 4"x4", 12 ply, 2 Sponges/Pouch, 25 Pouches/Box, 24 Boxes/Case, 1200 Sponges/Case	White	Sterile
D-4412	AMD-Ritmed Gauze Sponges, 4"x4", 12 Ply, 1 Sponges/Pouch, 100 Pouches/Box, 12 Boxes/Case, 1200 Sponges/Case	White	Sterile
D-4416	AMD-Ritmed Gauze Sponges, 4"x4", 16 Ply, 2 Sponges/Pouch, 25 Pouches/Box, 16 Boxes/Case, 800 Sponges/Case	White	Sterile
D-4800	AMD-Ritmed Gauze Sponges, 4"x4", 8 Ply, 2 Sponges/Pouch, 25 Pouches/Box, 24 Boxes/Case, 1200 Sponges/Case	White	Sterile

Folding diagram


8 Ply (C-Fold)




TECHNICAL DATA SHEET

Gauze Sponge Items # D-2120 to D-4800

12 Ply (C-Fold)



16 Ply (C-Fold)



Regulatory

FDA Class 1, FDA Product Code NAC
Health Canada Class 1

Storage Conditions

Store under ambient temperature and dry conditions

