

TECHNICAL DATA SHEET

Elastic Bandage Roll Items # A620 to A623

Intended Use

Elastic bandage used to provide flexible support and compression during activity. Ideal for strains and sprains. This type of product fits all parts of the body and can be washed and reused. To be used in health care facilities or over the counter.

Characteristic	Specification			
Description	Roll of woven elastic material with selvedge edges (self-finished) held by 2 metal clips. Extremities are fused.			
Material Composition	Spandex & Polyester. Latex-free.			
Color	Grey			
Dimension	Size	Item Code	Width (cm)	Length (m)
	2" x 5 YD	A620	5,0	4,5
	3" x 5 YD	A621	7,5	
	4" x 5 YD	A622	10,0	
	6" x 5 YD	A623	15,0	
Packaging	Wrap has a perforated line circling the roll in the middle of the length			
	Roll/Wrap	Wraps/Box	Boxes/Case	Rolls/Case
	1	10	5	50

Item Code	Description	Color	Specific Attributes
A620	AMD-Ritmed Elastic Bandage Roll, 2" x 5 YD (5,0 cm x 4,5 m), 1 Roll/Wrap, 10 Wraps/Box, 5 Boxes/Case, 50 Rolls/Case	Grey	Non-sterile, Individually Wrapped, Latex-free
A621	AMD-Ritmed Elastic Bandage Roll, 3" x 5YD (7,5 cm x 4,5 m), 1 Roll/Wrap, 10 Wraps/Box, 5 Boxes/Case, 50 Rolls/Case	Grey	Non-sterile, Individually Wrapped, Latex-free
A622	AMD-Ritmed Elastic Bandage Roll, 4" x 5YD (10,0 cm x 4,5 m), 1 Roll/Wrap, 10 Wraps/Box, 5 Boxes/Case, 50 Rolls/Case	Grey	Non-sterile, Individually Wrapped, Latex-free
A623	AMD-Ritmed Elastic Bandage Roll, 6" x 5 YD (15,0 cm x 4,5 m), 1 Roll/Wrap, 10 Wraps/Box, 5 Boxes/Case, 50 Rolls/Case	Grey	Non-sterile, Individually Wrapped, Latex-free

Visual





TECHNICAL DATA SHEET

Elastic Bandage Roll Items # A620 to A623

Regulatory

FDA Class 1, FDA Product Code FQM

Health Canada Class 1

Storage Conditions

Store under ambient temperature and dry conditions

Country of Origin

China